# BRAUN

## ThermoScan® 7

Ear thermometer

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IRT6520

English .....	2
Français .....	14
Español.....	25

## Contact Us

Questions or Comments? 

Call: 1-800-327-7226

Visit: [www.BraunHealthcare.com/contact-us](http://www.BraunHealthcare.com/contact-us)  
BraunHealthcare.com

Certain trademarks used under license from The Procter & Gamble Company or its affiliates. ThermoScan®, ExactTemp® and Age Precision™ are trademarks owned by Helen of Troy Limited.

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
IRT6520  
A001684R5  
15JUN22

## IRT 6520

- 1 Hygienic lens filter (Box of 20)
- 2 Probe tip
- 3 Probe
- 4 ExacTemp® Stability indicator light
- 5 Lens filter ejector
- 6 Display
- 7 Age Precision™ Technology and memory buttons
- 8 Power button
- 9 Start button
- 10 Battery door
- 11 Protective case

## Product description

The Braun ThermoScan® Ear thermometer has been carefully developed for accurate, safe, and fast temperature measurements in the ear.

The shape of the thermometer probe prevents it from being inserted too far into the ear canal, which can hurt the eardrum.

However, as with any thermometer, proper technique is critical to obtaining accurate temperatures. Therefore, read the instructions carefully and thoroughly.

## Indications for use / Intended use

The Braun ThermoScan® IRT6000 Series infrared ear thermometers are indicated for the intermittent measurement of human body temperature by consumers of all ages in a home use environment. The hygienic lens filter is used as a sanitary barrier between the thermometer and the ear canal.

## Warnings and Precautions

Please consult your physician if you see symptoms such as unexplained irritability, vomiting, diarrhea, dehydration, changes in appetite or activity, seizure, muscle pain, shivering, stiff neck, pain when urinating, etc. in spite of any color of the background light or absence of fever.

Please consult your physician if the thermometer shows fever (yellow or red background lights).

The operating ambient temperature range for this thermometer is 50 – 104 °F (10 – 40 °C).

**DO NOT** expose the thermometer to temperature extremes (below –13 °F / –25 °C or over 131 °F / 55 °C) or excessive humidity (15 - 95 % RH).

This thermometer must only be used with genuine Braun ThermoScan Hygienic lens filters (LF 40).

To avoid inaccurate measurements **ALWAYS** use this thermometer with a new, clean lens filter attached.

If the thermometer is accidentally used without a hygienic lens filter attached, clean the lens (see «Care and cleaning» section).

Keep hygienic lens filter out of reach of children.

This thermometer is intended for household use only.

Use of this thermometer is not intended as a substitute for consultation with your physician.

This thermometer is not intended for pre-term babies or small-for-gestational age newborn babies (birth to 4 weeks of age). Pre-term is defined as babies born before 37 weeks of gestation. Small-for-gestational-age is defined as a newborn baby, born at 37 weeks or later, with a weight below the 10th percentile for newborn babies of the same gestational age.

Nursing mothers should not take a temperature measurement on themselves or have their temperature measured by another individual while or immediately after nursing a baby, as nursing can affect body temperature.

This thermometer is not intended to interpret hypothermic temperatures. If the device displays a temperature of 97.5°F (36.4°C) or less with an LED backlight color of green or white, and the individual is exhibiting atypical symptoms or behaviors, contact your doctor or health care professional.

**DO NOT** allow children to take their temperature unattended.

**DO NOT** modify this equipment without authorization of the manufacturer.

People who are on antibiotics, analgesics, or antipyretics should not be assessed solely on temperature readings to determine the severity of their illness.

Temperature elevation may signal a serious illness, especially in neonates and infants, or in adults who are elderly, frail, or have a weakened immune system. Please seek professional advice immediately when there is a temperature elevation and if you are taking temperature on:

- Neonates and infants under 3 months (Consult your physician immediately if the temperature exceeds 99.4 °F or 37.4 °C).
- Individuals having diabetes mellitus or a weakened immune system (e.g., HIV positive, cancer chemotherapy, chronic immunosuppressant treatment, splenectomy).
- Individuals who are bedridden (e.g., paraplegia, quadriplegia, stroke, chronic illness, surgical recovery, nursing home patient, stroke, chronic illness, recovering from surgery).
- A transplant recipient (e.g., liver, heart, lung, kidney).

This thermometer contains small parts that can be swallowed or produce a choking hazard to children.

Keep out of reach of children.

## Body temperature


Normal body temperature is a range. It varies by site of measurement, and it tends to decrease with age. It also varies from person to person and fluctuates throughout the day. Therefore, it is important to determine normal temperature ranges. This is easily done using Braun ThermoScan. Practice taking temperatures on yourself and healthy family members to determine the normal temperature range.

Note: When consulting your physician, communicate that the ThermoScan temperature is a temperature measured in the ear and if possible, note the individual's normal ThermoScan temperature range as additional reference.

## Why measure in the ear?

The goal of thermometry is to measure core body temperature<sup>1</sup>, which is the temperature of the vital organs. Ear temperatures accurately reflect core body temperature<sup>2</sup> since the eardrum shares blood supply with the temperature control center in the brain,<sup>3</sup> the hypothalamus. Therefore, changes in body temperature are reflected sooner in the ear than at other sites. Axillary temperatures measure skin temperature and may not be a reliable indicator of core body temperature. Oral temperatures are influenced by drinking, eating and mouth breathing. Rectal temperatures often lag behind changes in core body temperature and there is a risk of cross-contamination.


## Age Precision™ Technology

Clinical research shows the definition of fever changes as newborns grow into children and children grow into adults.<sup>4</sup>

Select the appropriate age setting, take a temperature, and the display will glow green, yellow, or red to indicate whether the temperature is in a normal, fever, or high fever range, respectively. The device is only intended for the measurement of human body temperature. It is not intended to diagnose or monitor a serious medical condition or disease.

## How does Braun ThermoScan® work?

Braun ThermoScan measures the infrared heat generated by the eardrum and surrounding tissues. When placed in the ear, Braun ThermoScan continuously monitors the infrared energy and the ExactTemp® Stability Indicator displays a reading only if the thermometer detects that an accurate and stable temperature measurement has been taken.

<sup>1</sup> Guyton A C, Textbook of Medical Physiology, W.B. Saunders, Philadelphia, 1996, p 919.

<sup>2</sup> Guyton A C, Textbook of Medical Physiology, W.B. Saunders, Philadelphia, 1996, pp 754-5.

<sup>3</sup> Netter H F, Atlas of Human Anatomy, Novartis Medical Education, East Hanover, NJ, 1997, pp 63, 95.

<sup>4</sup> Herzog L, Phillips S. Addressing Concerns About Fever. Clinical Pediatrics. 2011; 50(#5): 383-390.

# How to use your Braun ThermoScan®

- 1** Remove the thermometer from the protective case.

- 2** Push the Power button .

During an internal self-check, the display shows all segments. Then the last temperature taken will be displayed for 5 seconds.


- 3** The hygienic lens filter indicator will blink to signal one is needed.

To achieve accurate measurements, make sure a new, clean, hygienic lens filter is in place before each measurement.

Attach a new hygienic lens filter by pushing the thermometer probe straight into the lens filter inside the box and then pulling out.


Note: The Braun ThermoScan will not work unless a hygienic lens filter is attached.


- 4** Select an age with the Age Precision™ Technology button  as shown. Press button to toggle through each age group.

**NOTE:** You must select an age in order to take a measurement.


**0-3 Months**


**3-36 Months**

**36 Months – adult**


The thermometer is ready to take a temperature when the display looks like the image below.

- 5 Fully insert the probe into the ear canal and point straight towards the opposite ear or temple.

Once positioned correctly, press and release the start button .

- 6 ExacTemp® Stability Indicator light will pulse while temperature measurement is in progress. The light will remain solid for 3 seconds to indicate that a successful temperature reading has been achieved.


**NOTE:** If the probe has been properly inserted into the ear canal during the measurement, a long beep will sound to signal the completed measurement.

If the probe has NOT been constantly placed in a stable position in the ear canal, a sequence of short beeps will sound, the ExacTemp® Stability Indicator light will go out and the display will show an error message (POS = position error).

See «Errors and troubleshooting» section for more information.

- 7 The confirmation beep indicates that an accurate temperature measurement has been taken. The result is shown on the display.


Age range	Green Normal temperature	Yellow Fever	Red High fever
 0*-3 Months	≥96.4 – ≤99.4 °F (≥35.8 – ≤37.4 °C)	X	>99.4 °F (>37.4 °C)
 3-36 Months	≥95.7 – ≤99.6 °F (≥35.4 – ≤37.6 °C)	>99.6 – ≤101.3 °F (>37.6 – ≤38.5 °C)	>101.3 °F (>38.5 °C)
 36 Months- adult	≥95.7 – ≤99.9 °F (≥35.4 – ≤37.7 °C)	>99.9 – ≤103.0 °F (>37.7 – ≤39.4 °C)	>103.0 °F (>39.4 °C)

\* 0 months is defined as having a gestational age of greater than or equal to 37 calendar weeks at birth.


- 8 For the next measurement, press Eject button  to remove, and discard used lens filter, and put on a new, clean lens filter.

**NOTE:** Defaults to last age setting used if you do not change.

The Braun ThermoScan Ear thermometer turns off automatically after 60 seconds of inactivity.

The thermometer can also be turned off by pressing and holding the Power button .

The display will briefly flash OFF and it will go blank.

## Temperature-taking hints

**ALWAYS** replace disposable lens filters to maintain accuracy and hygiene. The right ear measurement may differ from the measurement taken in the left ear. Therefore, **ALWAYS** take the temperature in the same ear. The ear must be free from obstructions or excess earwax build-up to take an accurate reading.

External factors may influence ear temperatures, including when an individual has:

Factor	Affects
Poor probe placement	✓
Used hygienic lens filter	✓
Dirty lens	✓

In the cases below, wait 20 minutes prior to taking a temperature.

Extreme hot or cold room temperatures	✓
Hearing aid	✓




Use the untreated ear if ear drops or other ear medications have been placed in the ear canal.

Use the ear that has been facing upwards if an individual has been lying down on their side with an ear face down on a pillow.

## Memory mode

- 1 This model stores the last 9 temperature measurements. To display the stored measurement, the thermometer must be turned on.
- 2 Press the **Mem** button. The display shows the temperature, and when releasing the **Mem** button, the stored temperature for that memory number is displayed, together with "M".
- 3 Each time the **Mem** button is pushed, the stored reading and an M is displayed to indicate each temperature measurement (e.g., M2).


Memory mode is automatically exited by not pressing the memory button for 5 seconds.

## Night light feature

The thermometer includes a convenient night light to illuminate the display in a dim environment.

The light will turn on when you press any button. It will stay on until the thermometer is inactive for up to 15 seconds even after temperature is taken.


## How to change temperature between °F and °C

US: Your Braun ThermoScan is shipped with the Fahrenheit (°F) temperature scale activated. If you wish to switch to Celsius (°C) and/or back from Fahrenheit to Celsius, proceed as follows:

CA: Your Braun ThermoScan is shipped with the Celsius (C) temperature scale activated. If you wish to switch to Fahrenheit (F) and/or back from Celsius to Fahrenheit, proceed as follows:

- 1 Make sure the thermometer is turned off.
- 2 Press and hold down the Power button . After about 3 seconds the display will show this sequence: °F / SET / °C / SET.
- 3 Release the Power  button when the desired temperature scale is shown. There will be a short beep to confirm the new setting, then the thermometer is turned off automatically.

## Care and cleaning

The probe tip is the most delicate part of the thermometer.

It must be clean and intact to ensure accurate readings.

If the thermometer is ever accidentally used without a hygienic lens filter, clean the probe tip as follows:

Very gently wipe the surface with a cotton swab or soft cloth moistened with alcohol. After the alcohol has completely dried out, you can put a new lens filter on and take a temperature measurement.

If the probe tip is damaged, please contact Consumer Relations.

Use a soft, dry cloth to clean the thermometer display and exterior. **DO NOT** use abrasive cleaners. **NEVER** submerge this thermometer in water or any other liquid. Store thermometer and lens filters in a dry location free from dust and contamination and away from direct sunlight.


Additional hygienic lens filter (LF 40) are available at most stores carrying Braun ThermoScan.

## Replacing the batteries

The thermometer is supplied with two 1.5 V type AA (LR06) batteries.

- 1 Insert new batteries when the battery symbol appears on the display.
- 2 Open the battery compartment. Remove the batteries and replace with new batteries, making sure the poles are in the right direction.
- 3 Snap battery door into place.




**Only discard empty batteries. To protect the environment, dispose of empty batteries at your retail store or at appropriate collection sites according to national or local regulations.**

**Do not mix alkaline, standard, or rechargeable batteries.**

## Calibration

The thermometer is initially calibrated at the time of manufacture. If this thermometer is used according to the use instructions, periodic re-adjustment is not required. If at any time you question the accuracy, please contact Consumer Relations.

## Errors and troubleshooting

Error message	Situation	Solution
	No hygienic lens filter is attached.	Attach new, clean lens filter.
	The thermometer probe was not positioned securely in the ear. An accurate measurement was not possible.	Take care that the positioning of the probe is correct and remains stable.
	Ambient temperature is not within the allowed operating range (50 – 104 °F or 10 – 40 °C).	Allow the thermometer to remain for 30 minutes in a room where the temperature is between 50 and 104 °F or 10 and 40 °C.
	Temperature taken is not within typical human temperature range (93.2 – 108 °F or 34 – 42.2 °C).	Make sure the probe tip and lens are clean and a new, clean lens filter is attached. Make sure the thermometer is properly inserted. Then, take a new temperature.
	HI = too high	
	LO = too low	
	System error – self-check display flashes continuously and will not be followed by the ready beep and the ready symbol or screen is blank.	Wait 1 minute until the thermometer turns off automatically, then turn on again.
	If error persists,	... reset the thermometer by removing the batteries and putting them back in.
	If error persists.	... replace the batteries with new batteries
	If error still persists,	... please contact the Consumer Relations.
	Battery is low, but thermometer will still operate correctly.	Insert new batteries.
	Battery is too low to take correct temperature measurement.	Insert new batteries.

## Product specifications

Displayed temperature range:	93.2 – 108 °F	(34 – 42.2 °C)
Operating ambient temperature range:	50 – 104 °F	(10 – 40 °C)
Storage temperature range:	-13 – 131 °F	(-25 – 55 °C)
Storage relative humidity:	Storage 15-95% RH (non-condensing)	
Operating relative humidity:	Operating 10-95%	
Display resolution:	0.1 °F or °C	

### Accuracy for displayed temperature range

35 – 42 °C (95 – 107.6 °F):	Maximum Laboratory Error	
Outside this range:	±0.4 °F	(±0.2 °C)
Clinical repeatability:	±0.5 °F	(±0.3 °C)
Reference Body Site:	±0.26 °F	(±0.14 °C)
Battery life:	Oral Equivalent	
Service life:	2 years in sleep mode / 350 measurements	
	5 years	

This thermometer is specified to operate at 1 atmosphere or at altitudes with an atmospheric pressure up to 1 atmosphere (700 - 1060 hPa).

Type BF Applied Part

Consult Instructions for Use

Operating temperature

Storage temperature

Keep dry

Attention, consult accompanying documents

15% Storage relative humidity

10% Operating relative humidity

Lot Number

Serial Number

Manufacturer

Date of Manufacture

Hygienic lens filter: Do not reuse

Subject to change without notice.

This appliance conforms to the following standards:

This infrared thermometer meets requirements established in ASTM Standard (E 1965-98) as it pertains to infrared thermometers. Full responsibility for the conformance of this product to the standard is assumed by Kaz USA, Inc. ASTM laboratory accuracy requirements in the display range of 98 to 102 °F (37 to 39 °C) for IR thermometers is  $\pm 0.4$  °F ( $\pm 0.2$  °C), whereas for mercury in-glass and electronic thermometers, the requirement per ASTM Standards E 667-86 and E 1112-86 is  $\pm 0.2$  °F ( $\pm 0.1$  °C).

EN ISO 80601-2-56: Medical electrical equipment – Part 2-56: Particular requirements for basic safety and essential performance of clinical thermometers for body temperature measurement.

EN IEC 60601-1: Medical electrical equipment – Part 1: General requirements for basic safety and essential performance.

EN ISO 14971: Medical devices – Application of risk management to medical devices.

EN ISO 10993-1: Biological evaluation of medical devices – Part 1: Evaluation and testing within a risk management process.

EN IEC 60601-1-2: Medical electrical equipment – Part 1-2: General requirements for basic safety and essential performance – Collateral standard: Electromagnetic compatibility – Requirements and tests.

EN ISO 15223-1: Medical devices. Symbols to be used with medical device labels, labelling and information to be supplied. General requirements.

EN IEC 60601-1-11: Medical electrical equipment – Part 1-11: General requirements for basic safety and essential performance – Collateral standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment.

This product conforms to the provisions of the EC directive 93/42/EEC.

MEDICAL ELECTRICAL EQUIPMENT needs special precautions regarding EMC. For detailed description of EMC requirements please contact your authorized local Service Center (See guarantee card).

Portable and mobile RF communications equipment can affect MEDICAL ELECTRICAL EQUIPMENT.

Please do not dispose of the product in the household waste at the end of its useful life.

To protect the environment, dispose of empty batteries at appropriate collection sites according to national or local regulations.

## Limited three year warranty

**You should first read all Instructions before attempting to use this product.**

- A. This 3 year limited warranty applies to repair or replacement of product found to be defective in material or workmanship. This warranty does not apply to damage resulting from commercial, abusive, unreasonable use or supplemental damage. Defects that are the result of normal wear and tear will not be considered manufacturing defects under this warranty. **KAZ IS NOT LIABLE FOR INCIDENTAL OR CONSEQUENTIAL DAMAGES OF ANY NATURE, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE ON THIS PRODUCT IS LIMITED IN DURATION TO THE DURATION OF THIS WARRANTY.** Some jurisdictions do not allow the exclusion or limitation of incidental or consequential damages or limitations on how long an implied warranty lasts, so the above limitations or exclusions may not apply to you. This warranty gives you specific legal rights, and you also may have other rights which vary from jurisdiction to jurisdiction. This warranty applies only to the original purchaser of this product from the original date of purchase.
- B. At its option, Kaz will repair or replace this product if it is found to be defective in material or workmanship.
- C. This warranty does not cover damage resulting from any unauthorized attempts to repair or from any use not in accordance with the instruction manual.

Call us toll-free at 1-800-327-7226 or visit: [BraunHealthcare.com/contact-us](http://BraunHealthcare.com/contact-us)

Please be sure to specify a model number.

**NOTE: IF YOU EXPERIENCE A PROBLEM, PLEASE CONTACT CONSUMER RELATIONS FIRST OR SEE YOUR WARRANTY. DO NOT RETURN THE PRODUCT TO THE ORIGINAL PLACE OF PURCHASE. DO NOT ATTEMPT TO OPEN THE THERMOMETER HOUSING YOURSELF, DOING SO MAY VOID YOUR WARRANTY AND CAUSE DAMAGE TO THE PRODUCT OR PERSONAL INJURY.**

**Guidance and manufacturer's declaration – electromagnetic immunity**

The IRT6520 is intended for use in the electromagnetic environment specified below.  
The customer or the user of the IRT6520 should ensure that it is used in such an environment.

Phenomenon	Basic EMC standard or test method	IMMUNITY TEST LEVELS
		HOME HEALTHCARE ENVIRONMENT
ELECTROSTATIC DISCHARGE	IEC 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air
Radiated RF EM fields <sup>a)</sup>	IEC 61000-4-3	10 V/m <sup>f)</sup> 80 MHz – 2,7 GHz <sup>b)</sup> 80 % AM at 1 kHz <sup>c)</sup>
Proximity fields from RF wireless communications equipment	IEC 61000-4-3	See RF wireless communication equipment immunity table below
RATED power frequency magnetic fields <sup>d) e)</sup>	IEC 61000-4-8	30 A/m <sup>g)</sup> 50 Hz or 60 Hz

a) The interface between the PATIENT physiological signal simulation, if used, and the ME EQUIPMENT or ME SYSTEM shall be located within 0,1 m of the vertical plane of the uniform field area in one orientation of the ME EQUIPMENT or ME SYSTEM.

b) ME EQUIPMENT and ME SYSTEMS that intentionally receive RF electromagnetic energy for the purpose of their operation shall be tested at the frequency of reception. Testing may be performed at other modulation frequencies identified by the RISK MANAGEMENT PROCESS. This test assesses the BASIC SAFETY and ESSENTIAL PERFORMANCE of an intentional receiver when an ambient signal is in the passband. It is understood that the receiver might not achieve normal reception during the test.

c) Testing may be performed at other modulation frequencies identified by the RISK MANAGEMENT PROCESS.

d) Applies only to ME EQUIPMENT and ME SYSTEMS with magnetically sensitive components or circuitry.

e) During the test, the ME EQUIPMENT or ME SYSTEM may be powered at any NOMINAL input voltage, but with the same frequency as the test signal (see Table 1).

f) Before modulation is applied.

g) This test level assumes a minimum distance between the ME EQUIPMENT or ME SYSTEM and sources of power frequency magnetic field of at least 15 cm. If the RISK ANALYSIS shows that the ME EQUIPMENT or ME SYSTEM will be used closer than 15 cm to sources of power frequency magnetic field, the IMMUNITY TEST LEVEL shall be adjusted as appropriate for the minimum expected distance.

**Guidance and manufacturer's declaration – electromagnetic emissions**

The IRT6520 equipment is intended for use in the electromagnetic environment specified below.  
The customer or the user of the IRT6520 should ensure that it is used in such an environment.

Emissions Test	Compliance	Electromagnetic environment – guidance
RF Emissions CISPR 11	Group 1	The ME equipment uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF Emissions CISPR 11	Class B	Complies
Harmonic emissions IEC 61000-3-2	Not Applicable	The ME equipment is solely battery powered.
Voltage fluctuations/ flicker emissions	Not Applicable	

Guidance and manufacturer's declaration – RF wireless communication equipment immunity						
Test frequency (MHz)	Band <sup>a)</sup> (MHz)	Service <sup>a)</sup>	Modulation <sup>b)</sup>	Maximum power (W)	Distance (m)	IMMUNITY TEST LEVEL (V/m)
385	380 – 390	TETRA 400	Pulse modulation <sup>b)</sup> 18 Hz	1,8	0,3	27
450	430 – 470	GMRS 460, FRS 460	FM <sup>c)</sup> ± 5 kHz deviation 1 kHz sine	2	0,3	28
710	704 – 787	LTE Band 13, 17	Pulse modulation <sup>b)</sup> 217 Hz	0,2	0,3	9
745						
780						
810	800 – 960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation <sup>b)</sup> 18 Hz	2	0,3	28
870						
930						
1 720	1 700 – 1 990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation <sup>b)</sup> 217 Hz	2	0,3	28
1 845						
1 970						
2 450	2 400 – 2 570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation <sup>b)</sup> 217 Hz	2	0,3	28
5 240	5 100 – 5 800	WLAN 802,11 a/n	Pulse modulation <sup>b)</sup> 217 Hz	0,2	0,3	9
5 500						
5 785						

<sup>a)</sup> For some services, only the uplink frequencies are included.

<sup>b)</sup> The carrier shall be modulated using a 50 % duty cycle square wave signal.

<sup>c)</sup> As an alternative to FM modulation, 50 % pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case.

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© 2022 All rights reserved.