# (•PainAwayPro•))

## with Heat

## Operating Manual

Item # RTLAGF-1000

Consult instructions for use

# Table of Contents

<u>Introduction</u> .....	1
<u>Indications for Use</u> .....	1
<u>Safety Warning</u> .....	2
<u>Contraindications</u> .....	2
<u>Warnings</u> .....	2
<u>Precautions</u> .....	3
<u>Adverse Reactions</u> .....	4
<u>Symbol and Title</u> .....	4, 5
<u>Environmental Condition for Transport and Storage</u> .....	6
<u>Electromagnetic Compatibility</u> .....	6, 7
<u>How the Device Works</u> .....	11
<u>Setup</u> .....	12
<u>Operating Instruction</u> .....	13
<u>Performance Specifications</u> .....	18
<u>Cleaning and Maintenance</u> .....	18
<u>Trouble Shooting</u> .....	18
<u>Recommended Use Positions</u> .....	19
<u>Warranty</u> .....	20
<u>Contact Information</u> .....	20

## **Introduction**

---

RTLAGF-1000 delivers electric pulses and heat generated to the user's body areas such as legs and feet through the electrodes. The portable and compact device has multiple modes of different pulse frequencies, covering Transcutaneous Electrical Nerve Stimulation (TENS) and Electrical Muscle Stimulation (EMS). It includes operating elements of ON/OFF button, intensity increase button, intensity decrease button, mode button, timer button, heat/temperature button, and could be attached and detached to electrodes through the magnetic connectors.

### **Indications for Use:**

#### **TENS (Modes 1, 2, 4, 5, 6, 8)**

To be used for temporary relief of pain associated with sore and aching muscles in the shoulders, waist, back, arms and legs, due to strain from exercise or normal household and work activities.

It is also intended for symptomatic relief and management of chronic, intractable pain and discomfort associated with arthritis.

#### **EMS (Modes 1, 3, 7)**

To stimulate healthy muscles to improve and facilitate muscle performance. To be used for the improvement of muscle tone and firmness and for strengthening muscles in the arms, abdomen, legs, and buttocks. Not intended for use in any therapy or for the treatment of any medical conditions or diseases.

It is also intended to temporarily increase local blood circulation in the healthy muscles of lower extremities.

### **Heating Mode**

Temporary relief of minor aches and pains.

## **Safety Warning**

---

Caution, consult accompanying documents

### **Contraindications:**

- ◆ Do not use this device on persons who have a cardiac pacemaker, implanted defibrillator, or other implanted metallic or electronic device because this may cause electric shock, burns, electrical interference or death.
- ◆ Do not use this device on persons whose pain syndromes are undiagnosed.

### **Warnings:**

- ◆ WARNING: Use carefully. May cause serious burns. Do not use over sensitive skin areas or in the presence of poor circulation. The unattended use of this device by children or incapacitated persons may be dangerous. To reduce the risk of burns, electric shock and fire this device must be used in accordance with the instructions.
- ◆ Do not crush the device or its electrode, and avoid sharp folds.
- ◆ Carefully examine the device and its electrodes, do not use if they show any sign of deterioration.
- ◆ Do not tamper with this device and its electrodes in any way. There are no user serviceable parts. If for any reason they do not function satisfactorily, return to the authorized service dealer where purchased from.
- ◆ The long-term effects of chronic electrical stimulation are unknown.
- ◆ Stimulation should not be applied over the carotid sinus nerves, particularly in persons with a known sensitivity to the carotid sinus reflex. Carotid sinus is located on both sides of the neck (*Refer to Page 20 for details*).
- ◆ Stimulation should not be applied over the front of neck or mouth. Severe spasm of the laryngeal and pharyngeal muscles may occur and the contractions may be strong enough to close the airway or cause difficulty in breathing. Stimulation over the neck can also cause adverse effects on heart rhythm or blood pressure.
- ◆ Do not apply stimulation across the patient's chest because the introduction of electrical current into the chest may cause rhythm disturbances to the patient's heart, which could be lethal.
- ◆ Since the effects of stimulation of the brain are unknown, stimulation should not be applied across the head and electrodes should not be placed on opposite sides of head.
- ◆ Stimulation should not be applied over swollen, infected, or inflamed areas or skin eruptions, e.g., phlebitis, thrombophlebitis, varicose veins, etc.

- ♦ Stimulation should not be applied over, or in proximity to, cancerous lesions.

### **Precautions:**

- ♦ Safety of stimulation use during pregnancy has not been established. Consult with your physician before use during pregnancy.
- ♦ Caution should be used for persons with suspected or diagnosed heart problems.
- ♦ Caution should be used for persons with suspected or diagnosed epilepsy.
- ♦ Caution should be used if you have any of the following:
  - If you have a tendency to bleed internally following an injury.
  - If you recently had surgery, or have ever had surgery on your back.
  - If areas of skin lack normal sensations, such as skin that is numb.
- ♦ Consult with your physician before use over the menstrual uterus.
- ♦ You may experience skin irritation or hypersensitivity due to the electrical stimulation or electrical conductive medium. The irritation can usually be reduced by using an alternate conductive medium, or alternate electrode placement.
- ♦ Do not use this device while driving, operating machinery, or during any activity in which involuntary muscle contractions may put the user at undue risk of injury.
- ♦ Keep this device out of reach of children.
- ♦ Do not use this device while sleeping.
- ♦ Do not use this device in high humidity areas such as, a bathroom.
- ♦ Stop using this device at once if you feel discomfort, dizziness or nausea, and consult your physician immediately.
- ♦ Do not attempt to move the electrodes while the device is operating.
- ♦ Do not apply stimulation of this device in the following conditions:
  - Across the chest because the introduction of electrical current into the chest may cause rhythm disturbances to the heart, which could be lethal.
  - Over open wounds or rashes, or over swollen, red, infected, or inflamed areas or skin eruptions (e.g., phlebitis, thrombophlebitis, varicose veins).
  - In the presence of electronic monitoring equipment (e.g., cardiac monitors, ECG alarms).

- On children or incapacitated persons.
- 
- ♦ Be aware of the following:
    - Consult with your physician before using this device.
    - This device is not effective for pain associated with Central Pain Syndromes, such as headaches.
      - This device is not a substitute for pain medications and other pain management therapies;
      - This device is a symptomatic treatment and, as such, suppresses the sensation of pain that would otherwise serve as a protective mechanism.
      - Stop using the device if the device does not provide pain relief.
      - Use this device only with the electrodes, and accessories recommended for use by the manufacturer.
    - ♦ Store the device away from high-temperature and direct-sunlight. Storage outside of stated storage temperature may result in measurement error or device malfunction.
    - ♦ Do not share the use of the electrodes with others; electrodes are intended for single person use.
    - ♦ Do not use the device while it's charging.
    - ♦ This device contains a lithium battery. If overheating of the device occurred during the charging, stop the charging or operation immediately and contact customer support.
    - ♦ Dispose of this battery-containing device according to the local, state, or federal laws.
    - ♦ Skin burns may occur, and check the skin under the electrode periodically.

## **Adverse Reactions**

- ♦ You may experience skin irritation and burns beneath the stimulation electrodes applied to the skin; check skin under electrodes frequently.
- ♦ You should stop using the device and should consult with your physicians if you experience adverse reactions from the device.


## Symbol and Title

Information essential for proper use shall be indicated by using the corresponding symbols. The following symbols may be seen on the device and labeling.

Symbol	Title
	Consult instructions for use
	Caution, consult accompanying documents
	Use by date
	Date of manufacture
	Batch code
	Catalog number
	Serial number
	Manufacturer
	Temperature limitation
	Humidity limitation
	Non-sterile
	MR Unsafe - keep away from magnetic resonance imaging (MRI) equipment
	Fragile, handle with care

	Keep away from rain
	Product packaging is able to be recycled

### Environmental Condition for Transport and Storage

 Fragile; handle with care

 Keep away from rain

 Product packaging is able to be recycled  
 Non-sterile

### Electromagnetic Compatibility

- 1) This product needs special precautions regarding electromagnetic compatibility (EMC) and needs to be installed and put into service according to the EMC information provided, and this unit can be affected by portable and mobile radio frequency (RF) communications equipment.
- 2) Do not use a mobile phone or other devices that emit electromagnetic fields, near the unit. This may result in incorrect operation of the unit.
- 3) Caution: This unit has been thoroughly tested and inspected to ensure proper performance and operation!
- 4) Caution: This machine should not be used adjacent to or stacked with other equipment and that if adjacent or stacked use is necessary, this machine should be

observed to verify normal operation in the configuration in which it will be used.

<b>Guidance and manufacture's declaration – electromagnetic emission</b>		
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should ensure that it is used in such an environment.		
<b>Emission test</b>	<b>Compliance</b>	<b>Electromagnetic environment – guidance</b>
RF emissions CISPR 11	Group 1	The device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emission CISPR 11	Class B	The device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

<b>Guidance and manufacture's declaration – electromagnetic immunity</b>			
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should ensure that it is used in such an environment.			
<b>Immunity test</b>	<b>IEC 60601 test level</b>	<b>Compliance level</b>	<b>Electromagnetic environment - guidance</b>
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floor are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2kV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	±1 kV differential mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% U <sub>T</sub> (>95% dip in U <sub>T</sub> ) for 0.5 cycle  40% U <sub>T</sub> (60% dip in U <sub>T</sub> ) for 5 cycles	<5% U <sub>T</sub> (>95% dip in U <sub>T</sub> ) for 0.5 cycle  40% U <sub>T</sub> (60% dip in U <sub>T</sub> ) for 5 cycles	Mains power quality should be that of a typical commercial or hospital environment. If the user of the device requires continued operation during power mains interruptions, it is recommended that the

	<p>70% <math>U_T</math> (30% dip in <math>U_T</math>) for 25 cycles</p> <p>&lt;5% <math>U_T</math> (&gt;95% dip in <math>U_T</math>) for 5 sec</p>	<p>70% <math>U_T</math> (30% dip in <math>U_T</math>) for 25 cycles</p> <p>&lt;5% <math>U_T</math> (&gt;95% dip in <math>U_T</math>) for 5 sec</p>	device be powered from an uninterruptible power supply or a battery.
Power frequency (50Hz/60Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE  $U_T$  is the a.c. mains voltage prior to application of the test level.

## Guidance and manufacturer's declaration – electromagnetic immunity

The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should ensure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3 V <sub>rms</sub> 150 kHz to 80 MHz	3 V <sub>rms</sub>	<p>Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p><b>Recommended separation distance</b></p> $d = 1,2\sqrt{P}$
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	$d = 1,2\sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = 2,3\sqrt{P} \quad 800 \text{ MHz to } 2,5 \text{ GHz}$ <p>Where <math>P</math> is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and <math>d</math> is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,<sup>a</sup> should be less than the compliance level in each frequency range.<sup>b</sup></p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the

device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the device.

b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

#### **Recommended separation distances between portable and mobile RF communications equipment and the device.**

The device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)		
	150 KHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 2.5 GHz
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

## How the Device Works

---


The device has multiple modes, covering TENS and EMS. If you are using the device for the first time, it is recommended that you start with the default Mode 1, which combines different pulse frequencies. Some modes are particularly effective for certain users, but you may need to select both the mode and the electrode position that are the best for you. You may need to try a few positions and modes before finding the one that suits you the best. The best choice of modes and positions may vary from one user to another, in spite of having the same type of symptom.

<b>Program</b>	<b>Pulse frequency (Hz)</b>
Mode 1: TENS and EMS	Combination of the following Mode 2-7
Mode 2: TENS	69
Mode 3: EMS	12.5-55.6
Mode 4: TENS	1.2
Mode 5: TENS	100 (On for 10 seconds and off for 2.5 seconds)
Mode 6: TENS	100 (On for 20 seconds and off for 1 seconds)
Mode 7: EMS	20
Mode 8: TENS	160

## Setup

---



Unpack the box of the product, and take the product and accessories out. The accessories may include electrodes, USB cable and AC adapter. The USB cable could be used to charge the control unit.

- ▶ Hold to turn on/off the device
- ▶ Press twice to lock/unlock the front buttons
- +
- 
- M Mode change
- T Timer change
- Heat +/- Turn on heat (Level 1), Level 2, and turn off heat
- ⎓ Charging port
- Connectors for stimulation and heat/temperature

## Operating Instruction

The following steps are used to guide the device operation, and the details about each step are listed in the following table.

<b>1<sup>st</sup> Step – Check the battery in the control unit</b>	<p>The control unit comes with a built-in rechargeable battery, and could be used as received. If the battery icon on the turned-on control unit keeps flashing, it means the battery is running out. Charge the control unit with the enclosed USB cable and AC adapter. USB cable could connect to any certified adapters. The battery icon flashes during charging, and becomes solid when the control unit is charged fully.</p>	 <p>Battery recharging</p>
<b>2<sup>nd</sup> Step - Install the electrode onto the control unit</b>	<p>Place the enclosed electrode onto the back side of the control unit using the magnetic snaps. Match up the silver magnetic snaps on the electrode pad with the silver magnetic snaps on the control unit, follow the same steps for the gold. This should be done before applying the device onto the skin of treatment areas.</p>	 <p>Electrode installment</p>
<b>3<sup>rd</sup> Step - Place the electrode-installed device on the stimulation-needed body area</b>	<p>Place the pad-installed device onto the treatment area (such as the shoulder). Press down firmly and ensure a firm contact with the body skin. <b>Note:</b> Keep the skin clean before placing the device.</p>	 <p>Place the device on the body area</p>
<b>4<sup>th</sup> Step - Press the “ON/OFF” to turn on the power</b>	<p>Holding the On/Off button could turn on the device, indicated by the lit-up display on the front of device. You could press and release the On/Off button to lock/unlock the buttons on the device front. When the buttons are locked, you could not operate them. However, you still could turn the device off by holding the On/Off button.</p>	
<b>5<sup>th</sup> Step - Select one of the stimulation modes</b>		

There are totally 8 preset stimulation modes. Press the “M” button to select a desired pulse mode. The mode selected will be shown on the display.

M		+
modo		
Time		
Intensity		
T	Heat/Temp	-

M Mode change

#### 6<sup>th</sup> Step - Choose the stimulation time

The default timer is 30 min, and pressing the “T” button could select a desired timer (30, 40, 50, 60, 10, and 20 min). The timer selected will be shown on the display.

M		+
modo		
Time		
Intensity		
T	Heat/Temp	-

T Timer change

#### 7<sup>th</sup> Step - Adjust the stimulation intensity

Pressing the “+” button could increase the intensity, while pressing the “-” button could decrease the intensity. The intensity change will be reflected on the display.

Note: With the increase of intensity (total of 20 levels), you may experience sensations like tingling, vibration, pain, etc. Therefore, gradually increase the intensity, and stop increasing when a comfortable level is reached.

#### 8<sup>th</sup> Step – Use heat/temperature on the device

Press the “Heat +/-” button for the first time to turn on the heat Level 1 of up to 104 °F, for the second time to adjust to Level 2 of up to 109.4 °F and for the third time to turn off the heat. The heat will be reflected on the display.

**Note:** The device is used for the electrical stimulation or heat; consult with your physician before using electrical stimulation and heat simultaneously.

#### 9<sup>th</sup> Step - Hold the “ON/OFF” again to turn off the power after done


When the countdown timer is up, the device will turn off automatically. The device could be also turned off by holding the On/Off button, indicated by the display light off.

M		+
modo		
Time		
Intensity		
T	Heat/Temp	-

+ Intensity increase  
- Intensity decrease

M		+
modo		
Time		
Intensity		
T	Heat/Temp	-

Heat +/- Heat on, increase, and off

On/Off

## Remote Operation

As shown above, the Electronic Pulse Stimulator could be operated independently by pressing its buttons (ON/OFF button, intensity increase/decrease buttons, mode/timer selection buttons and heat/temperature button). As the optional accessory, the remote control could provide a secondary operation way to wirelessly control the functions mentioned above. At any time, you could switch to manually operate the device itself by pressing its buttons to deal with any issues that may arise.

1. Follow the above instruction to install the electrode onto the device, place on the body area, and turn on the device.

2. Recharge the battery of the remote if needed by inserting the USB cable to the top of the remote, which could connect to the adapter included or any certified adapters.

3. Switch the ON/OFF button to turn on the remote, indicated by its first LCD display on.

**Note:** You could press and release the lock/unlock button on the remote to lock/unlock the buttons on the remote and device. When the buttons are locked, the lock icon on the remote appears, and you could not operate the buttons on the remote and device. However, you still could turn the remote and device off by switching the On/Off button.

4. Pressing and releasing the “+” button on the remote would increase the intensity of the device, while pressing and releasing the “-” button on the remote would decrease the intensity of the device. Each time when the intensity increase/decrease once, the intensity bar on the remote will accordingly increase/decrease once.

**Note:** With the increase of intensity (total of 20 levels), you may experience sensations like tingling, vibration, pain, etc. Therefore, gradually increase the intensity, and stop increasing when a comfortable level is reached.

5. Press and release the “HEAT” button on the remote for the first time to turn on the heat Level 1 of up to 104 °F, indicated by the first HEAT indicator flashing on the LCD display; press and release the “HEAT” button for the second time to adjust to Level 2 of up to 109.4 °F, indicated by the second HEAT indicator flashing on the LCD display; press and release the “HEAT” button for the third time to turn off the heat, indicated by the HEAT indicator stop flashing.

6. There are total of 8 preset stimulation modes. Press and release the “M/T” button on the remote to select a desired pulse mode, indicated by the flashing of the selected mode on the LCD display.
7. The default timer is 30 min. Pressing and holding the “M/T” button could select a desired timer in the order of 30, 40, 50, 60, 10, and 20 min, indicated on the LCD display.
8. When the countdown timer is up, the remote and device will turn off automatically. The remote and device could be also turned off by switching the On/Off button on the remote, indicated by the LCD display off.

**Note:** When the battery of the remote is running low, its battery icon will flash, meaning it is time to recharge the battery.

If using the device for the first time, you may start from the default Mode 1, which combines the different frequencies. If you expect to use the specific TENS or EMS mode, please refer to the following for details.

### **Use as TENS**

The Electronic Pulse Stimulator with Heat includes the following TENS frequencies, 69Hz (Mode 2), 1.2 Hz (Mode 4), 100Hz (Modes 5 and 6), and 160Hz (Modes 8). The device could be placed directly on the treatment site where you are experiencing pain (such as pain associated with sore and aching muscles and pain associated with arthritis). For example, if the treatment site is on the leg, place the TENS device above the pain site, and select one of the modes mentioned.

Arthritis involves inflammation on joints of the body, usually producing pain and stiffness. The electrode is placed on or near the area of the arthritis pain. The TENS mode of the device generates electrical pulses that are sent via the electrode for pain relief.

### **Use as EMS**

The Electronic Pulse Stimulator also includes the EMS frequency, 12.5-55.6Hz (Mode 3) and 20Hz (Mode 7). The device could be placed directly on the body area where you are expecting the muscle firming and strengthening. For example, if the body area is the arm, place the device on the arm site, and select Mode 3 or 7.

Also, the EMS mode could help temporarily increase local blood circulation in the healthy muscles of lower extremities through the electrode pad which is positioned on the lower extremity.

#### **Recommended Practice for Both TENS and EMS:**

- Duration suggested for each body area is 30 min and 2 times per day. Consult with your physician for longer and more frequent uses.
- Start from the lowest intensity and gradually adjust the intensity to a comfortable level at a scale from 1 to 20.
- Good skin care is important for a comfortable use of device. Be sure the treatment site is clean of dirt and body lotion.

- Keeping the electrode in the storage bag after use will extend its lifespan. The electrode pad is disposable and should be replaced when it loses the adhesiveness. To purchase additional electrodes, please contact the seller.

As shown in the above device operation, there is one kind of electrode pad used with the Pain Away Pro with Heat. The biocompatible electrode pad, consists of electrode gel, conductive film, heating layer, backing material, and connectors.

In addition, it is worth mentioning the button-lock feature is very practical, which could prevent any unintended operations during use. Pressing the On/Off button twice could activate/deactivate the button-lock feature. This feature would not affect the safety or effectiveness of the subject device, and the user could still turn off the device by holding the On/Off button. After selecting the appropriate mode, intensity, and timer, the user could activate the button-lock; when the timer is up, the device will turn off automatically.

Electrodes with a shelf life of 2 years are intended for single person use, and could be replaced. If the electrode no longer contacts well with your body skin, it is time to replace it. If needed, you could use a damp and clean cloth to wipe electrodes between uses. Please see the following for the Directions for Use, Removal and Storage of the electrode used for both TENS and EMS.

#### Directions for Use

1. Clean skin thoroughly prior to each application of electrodes, which will not contact well if any lotion, make-up, or dirt is left on the body skin.
2. Ensure the device is off before applying the electrode to it.
3. Apply the electrode-installed device firmly to the body skin.

#### Removal and Storage

1. Turn the device off before removing the electrode from the body skin.
2. Grab the edge of the electrode and remove from the body area.
3. When not in use and/or between each use, store the electrode in the re-sealable bag, out of direct sunlight.

#### Heating Mode

In addition to the TENS and EMS modes described above, the device has an additional mode of heat/temperature. When used in the heating function, the device will couple with an electronically controlled electrode pad or garment to provide automatic thermal heat to the skin with a range of 96.8 °F - 109.4 °F. This automatic thermal control mechanism could also avoid skin from thermal injury. The device is used for the electrical stimulation or heat; consult with your physician before using electrical stimulation and heat simultaneously.

## **Performance Specifications**

---

Power Source	3.7V Battery
Number of Output Modes	8 preset modes
Timer Range (minutes)	10-60
Dimensions (mm) [L x W x D]	89 x 77 x 18 mm
Waveform	Biphasic
Shape	Rectangular
Maximum Output Voltage	46V@500Ω
Maximum Output Current	92mA@500Ω
Maximum Pulse Duration	100µSec
Maximum Frequency	160Hz

## **Cleaning and Maintenance**

Please use a wipe to clean the device first, and then use a dry cloth to wipe it again. The electrodes coming with the device are disposable and should be replaced when they lose their adhesiveness. Contact the seller for replacements. Replacement Electrodes, Item# RTLAGF-1010, comes with 1 small and 1 large electrode.

## **Trouble Shooting**

---

If your device is not operating properly, please check below for common problems and suggested solutions. If the recommended action does not solve the problem, please contact the seller.

### **Stimulation is weak or non-existent**

- Be sure skin is clean and pads are firmly attached to skin.
- The battery is low and needs to be charged.

### **Device does not turn on**

- Check if battery is low and needs to be charged.

### **Electrodes are not sticky**

- The pads will lose adhesiveness over use. Please contact the distributor or dealer for replacements.

### **Skin turns red**

- Stop to treat another area.
- If problem persists, contact your physician.

### **Wireless control not connected**

- Restart the device and wireless controller to re-connect automatically.
- Switch to operate the device itself manually.

## Recommended Use Positions

---

SHOULDERS

BACK

BUTTOCKS

ABDOMEN

ARMS

LEGS

FEET

Never apply electrodes on the throat or both sides of the neck, where the carotid sinus nerves are located.

## **Drive PainAway Pro TENS with Heat Limited Warranty**

---

Please contact Drive DeVilbiss Healthcare in case of a claim under the warranty. If you must send in the unit, enclose a copy of your receipt and clearly state the defect. The following warranty terms apply:

1. The warranty period for the remote and receiver is two years from date of purchase. In case of a warranty claim, the date of purchase has to be proven by means of the sales receipt or invoice.
2. Replacement parts under warranty do not extend the warranty period for the device.
3. The following is excluded under the warranty:
  - a. All damage due to improper treatment. e.g. non-observance of the user instruction.
  - b. All damage which is due to repairs or tampering by the customer or unauthorized third parties.
  - c. Damage during transport from the manufacturer to the consumer or during transport to the service center.
  - d. The gel electrode pads and USB cable are subject to normal wear and tear.
4. Liability for direct or indirect consequential losses caused by the unit is excluded even if the damage to the unit is accepted as a warranty claim.

For more information, please contact the manufacturer.

---

### **Contact Information**

Drive DeVilbiss Healthcare  
99 Seaview Boulevard  
Port Washington, NY 11050

Phone: (516) 998-4600

Technical Support: (866) 621-4933

Consumer Customer Service: (866) 694-5085

Dealer and Provider Customer Service: (877) 224-0946

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