


COMFORTWIDE™ Series



Floor Level Height

BED LOWERS TO 3.6" OFF THE FLOOR AROUND
THE WIDTH OF A SMARTPHONE



AllCare BED®


- **Falls Risk** – Floor level (3.6" low height) with optional floor mat to address high risk resident falls
- **Bariatric Dignity** - Provide bariatric residents the same equipment (42" wide/600lb capacity) utilized throughout your facility
- **Custodial Care** – Enhance resident comfort (42" width) to take care of custodial residents.
- **Short Term Rehab** – Advanced features (OneButton Cardiac Chair positioning, trend/reverse trend), optional scale to manage higher acuity residents
- **Standardization** - Improve staff training as well as maintenance with ONE bedframe solution across your facility
- **Industry Leading Warranty** – Made in the USA with a best in class warranty (lifetime on the welds, 15 years on the frame, 5 years on the Linak motor/drive system)
- **NEW ENHANCEMENTS** - Clinical Contour reduces resident migration eliminating staff repositioning, backlit pendant with lockouts for easier/safer use, optional underbed light and optional staff controls, shorter head and footboards (AllcareC)

AllCare Bed SPECIFICATIONS

Deck Height Low	3.6" optional (7.25")
Deck Height High	20" optional (23.25")
Sleeping Surface	35"w, 39"w, 42"w x 80/84"
Maximum Safe Working Capacity	600 lbs
Product Weight	275 lbs



LIMITED WARRANTY

15 years on bed frame mechanics	1 year on wood, casters & handsets
5 years on drive system	Lifetime warranty on all welds



A SNF experienced a 100% Reduction in BedFalls when using the AllCare Frame and Safety Mat

ENDBOARD OPTIONS



ACCESSORIES



Scale



Pivoting Assist Bars



Staff Controls with lockouts



Trapeze Assembly



1/4 Length Side Rail



UnderBed Light



Safety Mat



Battery backup



LED Back-lit Pendant w/Lockouts

PROUDLY MADE IN THE U.S.A.

(877) 867-7365
MED-MIZER.COM



is at the forefront of innovating medical technology for post acute care. We are focused on partnering with you to help address healthcare challenges in the following key areas—




FALLS PREVENTION



BARIATRIC CARE



MOBILITY



REHABILITATION



WOUND CARE