

COMFORTWIDE™ Series



Floor Level Height

BED LOWERS TO 3.6" OFF THE FLOOR AROUND
THE WIDTH OF A SMARTPHONE



Shown with optional
Shaker Style Boards

AllCare BED®

- **Falls Risk** – Floor level (3.6" low height) with optional floor mat and advanced falls notification system to address high risk resident falls
- **Bariatric Dignity** - Provide bariatric residents the same equipment (42" wide/600lb capacity) utilized throughout your facility.
- **Custodial Care** – Enhance resident comfort (42" width) and height (up to 30") to take care of custodial residents.
- **Short Term Rehab** – Advanced features (cardiac chair positioning, trend/reverse trend), optional scale, optional bed exit to manage higher acuity residents.
- **Standardization** - Improve staff training as well as maintenance with ONE bedframe solution across your facility.
- **Industry Leading Warranty** – Made in the USA with a best in class warranty (lifetime on the welds, 15 years on the frame, 5 years on the Linak motor/drive system)

AllCare Bed SPECIFICATIONS

Deck Height Low	3.6" optional (8.5")
Deck Height High	25" optional (30")
Sleeping Surface	35"w, 39"w, 42"w x 80/84"l
Maximum Safe Working Capacity	600 lbs
Product Weight	275 lbs




LIMITED WARRANTY

15 years on bed frame mechanics	1 year on wood, casters & handsets
5 years on drive system	Lifetime warranty on all welds

A SNF experienced a 100% Reduction in BedFalls when using the AllCare Frame and Safety Mat

ENDBOARD OPTIONS



ACCESSORIES

Scale

Pivoting Assist Bars

Advanced Falls Notification System

Trapeze Assembly

1/4 Length Side Rail

Safety Mat

Battery backup

PROUDLY MADE IN THE U.S.A.

(877) 867-7365
MED-MIZER.COM



is at the forefront of innovating medical technology for post acute care. We are focused on partnering with you to help address healthcare challenges in the following key areas—




FALLS PREVENTION



BARIATRIC CARE



MOBILITY



REHABILITATION



WOUND CARE