


THE MIGHTY SLING OPTION CAN BE USED ON  
EITHER THE MIGHTY 400 OR THE MIGHTY 600  
AS WELL AS THE PORTABLE MIGHTY VOYAGER.

SLING HEIGHT AND DEPTH DIMENSIONS ARE  
APPROXIMATE AS THE SLING CAN BE ADJUSTED  
TO FIT THE USERS NEEDS.



**TOLERANCES**

DIMENSIONS ARE IN INCHES  
UNLESS OTHERWISE SHOWN  
FRACTIONAL: .....  $\pm 1/16$   
ANGULAR: .....  $\pm 1^\circ$   
TWO PLACE DECIMAL: ....  $\pm .01$   
THREE PLACE DECIMAL: ..  $\pm .005$

**SLING ATTACHMENT - MIGHTY**  
**U.S. PATENT NO. D507,769S**  
**PART NO. F-MTYSLS**

**MATERIAL: 304 STAINLESS STEEL**  
**FINISH: POWDER COAT - WHITE**  
**APPROX. WEIGHT: 187.30 LBS.**



**9889 GARRYMORE LANE**  
**MISSOULA, MT 59808**  
**(406) 549-0769**  
**FAX (406) 549-2602**  
**www.AquaCreek.com**

**PROPRIETY AND CONFIDENTIAL INFORMATION**

THIS DRAWING AND ALL INFORMATION CONTAINED HEREIN  
IS THE SOLE PROPERTY OF AQUA CREEK PRODUCTS LLC. ANY  
REPRODUCTION IN PART OR AS A WHOLE WITHOUT WRITTEN  
PERMISSION FROM AQUA CREEK IS PROHIBITED.

\*SPECIFICATIONS ARE SUBJECT TO CHANGE WITHOUT NOTICE

REVISIONS			
REV.	DESCRIPTION	DATE	APPROVED
P1	PRODUCTION RELEASE	4/9/2021	J.MAXWELL