



ReOrder #	Description	Size	Units
909368	Auxiliary Connector Cable	3 ft Long	5pk

Purpose

The Auxiliary Cable is required to interface the following sensors with any of the ProSeries alarms.

Skil-Care Sensor Chart - Table 1

ReOrder #	Description	Size	Units
909391	UnderSeat Sensor	Fits 18"-22" Wheelchairs	Each
909392	Seat Belt w/Buckle Sensor	50"L	Each
909394	Seat Belt w/Buckle Sensor	70"L	Each
909530	Gel-Foam Cushion Sensor	16" x 16" x 2"	Each
909531	Gel-Foam Cushion Sensor	18" x 16" x 2"	Each
909532	Gel-Foam Cushion Sensor	20" x 16" x 2"	Each
909535	Leveling Pad Sensor	Fits 16"-18" Wheelchairs	Each

This cable is included with any of the chair alarm systems or ChairPro alarm (model# 909367) or the MultiPro alarm (model# 909510). In the event if any of the above sensors (listed in table 1) are used with the BedPro alarm, this cable must be purchased separately (model# 909368). This cable isn't required for any of the flat chair or mattress sensor pads or the Lift-Off Lap Cushion.

Caution

If you have purchased a ChairPro or MultiPro alarm do not discard cord since it may be needed at a later time.