

INVACARE® OVERBED TABLES



Invacare® Overbed Tables are practical, easy-to-use units suitable for both home and institution. The table height easily adjusts for comfort while reading, writing or eating. These tables feature an attractive wood-grained laminate top in your choice of finishes with a T-molded edge that is designed to resist damage. The heavy-gauge tubular welded steel construction is designed to ensure durability and the casters ensure smooth mobility of the table in any direction.



Invacare® Overbed Table

Model # IH5135

With a travel range of 30" to 45", this heavy duty overbed table features durable, robust construction with attractive, functional design. The extra-wide base is designed to accommodate most wheelchairs and resident seating to provide a supportive surface for eating, reading or most any activity.

The table top easily raises with the slightest upward touch and locks securely into place. A simple squeeze of the release lever allows the table top to lower to just the right height. Premium high pressure laminate tops provide improved durability and are available in your choice of eight finishes.

Large 2" casters are standard with this table and provide improved maneuverability.



Invacare® Overbed Table

Model # IH1010

Invacare's Overbed Table is an economical choice with advanced features. Ideal for home and institutional settings, these overbed tables provide a chrome U-base design that allows for positioning around most wheelchairs.

The table top easily raises from 28" to 45" with the slightest upward touch and a simple squeeze of the release lever allows the table top to lower. Premium high pressure laminate tops provide improved durability and are available in your choice of eight finishes.

Standard 1-1/2" casters allow for easy relocation.

Invacare Corporation
www.invacare.com

USA
One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada
570 Matheson Blvd. E.,
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

©2016
Invacare Corporation.
All rights reserved.
Trademarks are
identified by the
symbols ™ and ©.
All trademarks are
owned by or licensed
to Invacare Corporation
unless otherwise
noted. Specifications
are subject to change
without notification.
Form No. 16-065
160322

Available Finishes



Model	Dimensions Table Top	Dimensions Base	Height Adjustment	Frame Construction	Weight	Weight Capacity	Limited Warranty
IH5135	32" L 18" W .75" D	28" L 15.5" W	30"- 45"	11 & 12 Gauge	54 lb.	100 lb.	1 Year
IH1010	30" L 15" W .75" D	24" L 15.5" W	28"- 45"	16 & 17 Gauge	19 lb.	50 lb.	1 Year

