



Completely Coordinated.  
Completely Cost Effective.  
Completely Easy.

**RAINFOREST FOLLIES SCALE COMPLETE**



# Imagination Series, Rainforest Follies, Pediatric Scale Table

## MODEL 7832

### TABLE FEATURES

- Clinton original illustrations
- Easy to clean laminate surfaces
- All graphics are permanently embedded in the high pressure laminate and will not peel off
- Graphics are on one long side and two short ends. Back is slate gray. Only manufactured as pictured.
- Comfortable work height of 36"
- Top is rimmed on three sides for safety
- Top's front and back rim edges are beveled and the corners rounded for safety
- Easy-clean, replaceable center pad
- Built-in infantometer with folding ends for measuring height
- Easy access to weight, measurement and storage
- Choose to access concealed paper dispenser from either end
- Recessed front toe kick
- Easy-clean all laminate base
- Ample 2 door and 1 drawer storage
- Soft-close, door hinges



### MODEL 7832 SPECIFICATIONS

43"L x 24"W x 38.5"H

### SCALE FEATURES

- Dependable, accurate Doran® scale
- Easy-read digital display
- Auto calibration on start-up
- Auto shut-off
- Electric or battery operation (6 AA batteries required)
- Weights in increments of lbs., ozs. or kilos
- Minimum registered weight is 5 lbs.—Maximum weight for scale is 40 lbs.

California residents please be advised of the following, pursuant to Proposition 65:

**WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

1

MANUFACTURED IN THE  
**USA**

### Options



78-P Parrot Scale Pal



066 Door & Inside Latch Combo




055 Drawer Locks

### Upholstery

Gunmetal

## Features



### Infantometer

Built-in infantometer with folding ends for measuring height



### Paper Dispenser

Choose to access concealed paper dispenser from either end



### Rimmed Top with Beveled Edges

Top is rimmed on three sides for safety and the front and back rim edges are beveled and the corners rounded for safety



### Padded Top

Easy-clean, padded, upholstered top for comfort



### Scale Head

Easy-read digital display



## Rainforest Follies Cabinets

### MODEL 6132-BW

#### "Rainforest Follies" Base Cabinets

4 doors with adjustable shelves and 3 drawers.  
Shipped in one unit.

#### "Rainforest Follies" Wall Cabinets

4 doors with adjustable shelves. Shipped in one unit.

Base & wall cabinets are only sold in sets

- Clinton original illustrations
- Pre-assembled units save time
- Easy-clean all laminate construction
- All graphics are permanently embedded in the high pressure laminate and will not peel off
- Tough, chip-resistant PVC edges
- Assembled with steel fasteners
- 4 adjustable leveling feet on each unit
- 3-way adjustable, soft-close 110° hinges
- Flat panel drawer fronts with metal sides that hold up to 75 lbs.
- Recessed toe kick
- Easy hook and rail installation for upper cabinets
- Slate Gray edgebanding on all edges
- Customizable with 4 options
- Easily moves for relocation
- Drawers can be added to any base unit
- Doors and drawer can be moved to accommodate sink location
- Shown with optional sink and faucet



#### SPECIFICATIONS

BASE: 60"L x 18"D x 35"H

WALL: 60"L x 12"D x 24"H

California residents please be advised of the following, pursuant to Proposition 65:

**WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

1



#### Options



066 Door & Inside Latch Combo



055 Door Locks




044 Add a Cabinet Drawer



22 Sink & Faucet

# 10' X 10' FLOOR PLAN

*Optional Items in Red*



1 Square = 1 foot

**NOTE:** These are suggested floor plans. Please check local building and fire codes before beginning any construction or renovations.




*Walls*  
2046-20 Garden Green  
561 Pistachio



*Flooring*  
Rejuvenations™ Timberline®  
Maple 37350 you better

**OPTIONAL ACCESSORIES** *(Shown in photo on front)*

- **03-CC** Tree Wall Sticker
- **T-70** Clear Labeled Sundry Jars (5 in set)
- **GS-3020** Double Stainless Steel Glove Box Holder
- **C-40** Side Chairs in Royal Blue upholstery (3RB)



800-441-9131 • 717-848-3519  
Toll-Free Fax 866-522-1233  
Fax 717-843-5871  
[www.clinton-ind.com](http://www.clinton-ind.com)