

AeroCarbon™ Ultra-Lite Rollator

Premium Performance Without the Premium Price



Bring lightweight carbon fiber into your everyday lineup with the affordable **AeroCarbon Ultra-Lite Rollator**! Its easy-to-maneuver frame and full suite of premium features, combined with MAP pricing, deliver strong margins for you and standout value for your customers.

Effortless Maneuverability, Travel & Storage

Easy to lift, carry, and maneuver at just 11.5 pounds. Stores compactly in closets, trunks, and more thanks to a slim fold.

Durable & Easy to Maintain

Provides structural integrity without the bulk with carbon fiber construction. Reduces common fracture points with a weld-free frame and supports simplified maintenance with fewer replacement parts.

Style Customers Love

Converts even hesitant mobility aid users with a sleek design and sophisticated matte finish.

Confidence & Comfort

Offers stable performance indoors and outdoors with 9.5" front wheels and 8" rear wheels, plus snag-free use with enclosed brake cables. Customize support with adjustable handle height and 17.5" wide seat with ergonomic backrest.

Product Specs

Frame Material	Carbon Fiber
Frame Color	Gray, Black
Wheel Size	9.5" front 8" back
Handle Height	31.7" – 37.6"
Product Dimensions	23.6" (W) x 26.8" (D)
Seat Dimensions	17.5" (W) x 9" (D)
Seat to Floor Height	23"
Product Weight	11.5 lb.
Weight Capacity	330 lb.
Shipping Weight	15.3 lb.
HCPCS Code	E0143+E0156
Warranty	Limited Lifetime

Item #	Description	UPC	UOM
10278GY	AeroCarbon™ Ultra-Lite Rollator, Gray	822383030883	1/ea
10278BK	AeroCarbon™ Ultra-Lite Rollator, Black	822383031194	1/ea