

## limited two year warranty

Your Drive branded product is warrantied to be free of defects in materials and workmanship for two years from the date of purchase for the original consumer purchaser.

This device was built to exacting standards and carefully inspected prior to shipment. This Limited Two Year Warranty is an expression of our confidence in the materials and workmanship of our products and our assurance to the consumer of years of dependable service.

This warranty does not cover device failure due to owner misuse or negligence, or normal wear and tear. The warranty does not extend to non-durable components, such as rubber accessories, casters, and grips, which are subject to normal wear and need periodic replacement.

If you have a question about your Drive device or this warranty, please contact an authorized Drive dealer.



© 2014 Medical Depot, Inc. All rights reserved.

Drive is a trademark of Medical Depot, Inc.

Port Washington N.Y. 11050 USA Made in China

## legs with casters

item #s FC 2029N & FC 4029N

for use with first class chair  
(item # FC 2000N & FC 4000N)

WenzeliteRe/hab<sup>TM</sup>  
drive



EU Authorized Representative  
Drive Medical LTD  
Ainley's Industrial Estate  
Elland, West Yorkshire,  
United Kingdom HX5 9JP

[www.drivemedical.com](http://www.drivemedical.com)

## legs with casters



Item # FC 4029N Shown

## assembly instructions

1. Loosen knobs (A) and remove all four legs.
2. Insert legs with anti-tipper in the back and the other two in the front. Adjust height as needed.



**Wenzel***ite*Re/hab  
**drive**

**Wenzel***ite*Re/hab  
**drive**