

# TRUE BARIATRIC PATIENT LIFT WITH DUAL MOTORS



## FREE SP:RIT

### 850 PRO

TRUE BARIATRIC 850 LB. PATIENT LIFT



- ULTRA POWERFUL
- HEAVY DUTY
- TRUE BARIATRIC
- SECURE AND STURDY

PROVIDES A SAFE ALTERNATIVE TO  
MANUAL LIFTING FOR BOTH THE  
CAREGIVER AND PATIENT

**850 Lbs.**  
Weight  
Capacity



#### ● DUAL MOTORS

Raises/Lowers the patient from floor to 77" high

#### ● POWER BASE

Expands from 26"-50" ultra wide

**WEIGHT CAPACITY:** 850 Lbs.

**HEIGHT AT SLING HOOK UP:** 31" - 77" high

**BASE WIDTH:** Open- 50" / Closed- 26"

**BED CLEARANCE:** At least 4.5"

**CASTER SIZE:** Rear- 4" locking / Front- 3" swivel

MODEL # **FS850**



# TRUE BARIATRIC PATIENT LIFT WITH DUAL MOTORS



## FREE SP:RIT

### 850 PRO

TRUE BARIATRIC 850 LB. PATIENT LIFT



- Dual motors power the lift range and legs to open and close
- Ultra Wide leg opening easily wraps around bariatric wheelchairs and equipment
- 50" extra long leg length for increased stability
- Oversized EZ Grip ergonomic handle for increased control when moving the lift
- Emergency stop button locks the lift in place
- CPR function lowers and raises the lift when the battery is discharged
- Total pinch point protection
- Hanging bar rotates 360 degrees in all directions for easy patient handling



Large oversized hooks with safety clip hold the sling securely in place



Easy Roll Dual Front Castors



**BARIATRIC DIGITAL SCALE (OPTIONAL)**

Large easy to read display



**LINAK**

**LINAK BATTERY AND MOTORS**

- Quick release battery eliminates downtime
- Battery can be replaced immediately with separate charging base
- Battery can also be charged on board
- LCD Display indicates battery charge level

