

WORLD'S LIGHTEST SCOOTER!

NEW

Feather™ SCOOTER 27X

At just 27 lbs.*, the Feather™ Scooter 27X is the lightest mobility scooter on the planet.

Easily take the NEW Feather™ Scooter 27X with you everywhere! Weighing only 27 pounds* and folding to a compact 16" x 24" - this scooter is the perfect travel companion. Its compact size also makes for easy storage when not in use.



FEATURES:

- only 27 lbs*
- 6 mile drive range
- one-fold frame
- airline compliant battery^

SPECIFICATIONS

Weight Capacity	250 lbs.	Ground Clearance	2"
Product Weight w/out batteries	27 lbs.	Front Wheel Size	6"
Product Weight w/batteries	30 lbs.	Rear Wheel Size	6"
Seat Width	16"	Turning Radius	51"
Seat Depth	11"	Max Climbing Angle	6 degrees
Seat-to-Floor Height	19"	Max Speed	4 mph
Overall Width	16"	Driving Range	6 Miles
Overall Length	24"	Battery Type	25.2V 6.6 AH
Overall Height	32"	Warranty on Frame	Three Year Limited
Folded Depth	16"	Warranty on Wearable Parts & Electronics	One Year Limited
Folded Height	10"	Warranty on Batteries	Six Months Limited

Lightweight



Foldable



Locking Brakes



Arrives Assembled



*weight without battery. ^removable lithium battery is airline compliant. Check with individual carrier for guidelines.