

BARIATRIC GERIATRIC RECLINER CHAIR



**WEIGHT CAPACITY:
500 LBS.**

**3 POSITION - UPRIGHT, DEEP
RECLINE, ELEVATED LEG/FOOT**

- Built-in Headrest
- Extra wide seat width: 25"
- Weight capacity: 500 lbs.
- Durable 5" castors (2 locking)
- Large fold out tray included

MODEL # **GC500**



BARIATRIC ROLLATOR



**WEIGHT CAPACITY:
500 LBS.**

- Soft padded oversized seat
- Weight capacity: 500 lbs.
- Castors: Durable 8" with soft-grip tires
- Basket included

MODEL # **RB500**