



With a high capacity and mechanical height rod, this scale properly weighs and measures the heaviest ambulatory patients. The 350° swivel head makes it easy for clinicians to view measurements and interact with the scale.

Advantages

- EMR Capable**
Scale interfaces with a computer, EMR software or other electronic device via Health o meter® Professional Connectivity Solutions
- Value-Added Features**
Operating functions (LB/KG Conversion, LB/KG Lock Out, BMI, Zero, Tare, Hold/Release, Auto Zero, Auto Off) included at no additional cost
- 350° Swivel Display**
Eye-level rotating display can be read from almost any angle and allows for ease of use and patient privacy
- Calculates Body Mass Index (BMI)**
Scale functions include the ability to calculate the patient's BMI
- Built-in Height Rod**
The easy-to-use height rod provides height measurements from 23 5/8" to 84" /60cm – 213cm



350° swivel head display



Specifications

- Capacity: 600 lb / 272 kg
- Resolution: 0.2 lb / 0.1 kg
- Platform Size (w x h x d):
14 1/4" x 14 1/4" x 2 5/8"
362 mm x 362 mm x 67 mm
- Product Footprint (w x h x d):
14 3/4" x 23 1/4" x 54 1/2"
375 mm x 591 mm x 1384 mm
- Product Weight: 24 lb / 11 kg
- Height Rod: 23 5/8" - 84" / 60 cm - 213 cm
- Display: 1" / 25 mm LCD, 350° Swivel
- Wheels: 2
- Connectivity: USB, optional Pelstar wireless technology
- Power Source: 100-240V adapter (ADPT50) included,
6 C-cell batteries (not included)
- Functions: LB / KG Conversion, LB / KG Lock Out , BMI, Zero,
Tare, Hold / Release, Auto Zero, Auto Off
- 2 year limited warranty

Model	Description
597KG	Heavy duty eye level digital scale, kilograms only
SS-597KL	ScaleSurance 2 year extended warranty for 597 models*
ADPT50	Adapter included, item # for replacement or spare
PLUGSET50	Plug adapters for use outside North America
PELSTARONE-P	Wireless thermal printer

* Not available in all countries