

or more than 20 years, Invacare has set the standard in manual wheelchairs for the HME industry.

But our customers asked us to raise that standard to even greater heights. We listened and said, "Yes, we can build an even better manual wheelchair."

As a result, Invacare has developed a totally new line of manual wheelchairs – the Invacare® IVC™ Manual Wheelchair Series, focusing on Interchangeability, Value and Complexity Reduction:

- IVC Manual Wheelchairs are Interchangeable and compatible throughout the entire current line of wheelchairs. In addition, these chairs offer interchangeability with the Tracer® and 9000 series too.
- The IVC line provides increased Value, with more features, more options, and more functionality, all at a better price.
- At the same time, providers asked us to make the line simple and easy to understand. We listened and responded. The new IVC wheelchairs represent a conscious effort at Complexity Reduction consolidating 18 different frames down to 10 frames, covering seven different Medicare codes. Fewer models make it easier to make the right choice for each consumer without sacrificing the features and benefits that make the chair a perfect fit for each consumer.

In addition, we have responded to your feedback with a new "no tools adjustment" front-rigging design, a new permanent arm design, a reclining-back kit and much more.

The new disk wheel lock is designed to provide longer wear and greater durability. This new patent ed design will reduce maintenence and service repair associated with the old bolted assembly.

Invacare has expended considerable resources in researching and designing the IVC line; however, what we didn't do was check with all those off-brand, imported chairs to see if our new line will be interchangeable with them. So why take a chance? Stick with the industry leader: Invacare.



Bariatric/Heavy-Duty Designation

Invacare products with the Bariatric/Heavy-Duty icon are designed for your bariatric patients whose height and weight require the use of a heavy-duty model. These products provide the same quality that all Invacare products possess. They also demonstrate our commitment to this product category.

## **Table of Contents**

Tracer® EX2	3
Transport	4
Tracer® SX5	5
Tracer® SX5 Full Recliner	6
9000 SL	7
9000 XT	8 – 9
9000 XT Full Recliner	10
9000 Topaz™	П
9000 XDT	12
Tracer® IV	13
Invacare® 900 Series	14
Invacare® 900 Bariatric Recliner	15
9000 Jymni™	16
9000 Jymni™ Recliner	17
Ramps	18
Accessories	19 – 22
Color Guide	23
Accessory Availability Guide	24 – 25
Dimensions Chart	26 – 27





Tracer offset permanent arms allow for full seat width between arms.

## Invacare® IVC™ Tracer® EX2 Wheelchair

The Invacare IVC Tracer EX2 wheelchair sets the new standard in manual wheelchairs. The Tracer EX2 combines the design and technology of the Invacare 9000 series and incorporates it into the Invacare Tracer series. A true dual axle position allows for repositioning the 24" rear wheels and 8" casters for adult and hemi seat-to-floor heights. Offset permanent arms offer full seat width up to 20" with a permanent arm, hemi spaced footrest, full range of accessories and a limited lifetime warranty. The new design also makes it possible to interchange components with the 9000 series chairs.

#### **FEATURES**

- Durable, low-maintenance, triple chromeplated, carbon steel frame is long-lasting
- Dual axle positions allow variation of seatto-floor heights
- Offset permanent arms available in 16", 18" and 20" widths
- 14-gauge crossbraces add strength and durability
- Urethane rear tires mounted on "no flex" wheels offer superior performance
- Low-maintenance, dual-sealed precision bearings
- Heavy-duty inner liners keep seat and back from stretching
- Flame-retardant upholstery resists mildew and bacteria, as well as cold-crack and UV damage
- · Upholstery is dual-embossed
- Meets RESNA standards for manual wheelchairs in product performance

Dimensions	
Seat sizes	
Widths	16", 18", 20"
Depth	16" standard
Seat-to-floor height	17.5", 19.5"
•	15.5", 16.5" 18.5"
	with different
	wheel and caster packages
Back height	16" fixed
Front riggings	
Hemi swingaway footrest	17"-21"
Hemi swingaway padded elevating legrest	17.5"-21"
Articulating elevating legrest	17.5"-21"
(Footplates available in aluminum or composite)	
Arm styles	
Permanent full-length, padded	
Full-length fixed height, padded, removable	
Desk-length fixed height, padded, removable	
(All arms are conventional)	
Weight capacity	250 lb.
Product weight (TREX2/WD86/28)	36 lb.
Limited warranty	lifetime on frame
·	and crossbraces
	3



# Invacare® IVC™ Transport Wheelchair

The Invacare IVC Transport wheelchair is designed for safe and efficient patient transport. The Transport has the added advantage of sharing interchangeable parts with the other wheelchairs in the Tracer series. The Transport is compact and lightweight.

#### **FEATURES**

- Durable, low-maintenance, triple chromeplated, carbon steel frame is long-lasting
- · Fold-down back makes storage convenient
- Welded caster forks supply durability and long, maintenance-free life
- Heavy-duty inner liners keep seat and back from stretching
- Flame-retardant upholstery resists mildew and bacteria, as well as cold-crack and UV damage
- Padded armrests provide additional patient comfort

Dimensions	
Seat sizes	
Widths	16", 18"
Depth	16"
Seat-to-floor height	21"
Back height	19" fold-down
Front riggings	
Swingaway footrest	17"-21"
Swingaway elevating legrest	17.5"-21"
Arm styles	
Permanent full-length, padded	
(All arms are conventional)	
Weight capacity	250 lb.
Product weight (TRAN86)	27 lb.
Limited warranty	3 years on frame and crossbraces



## Invacare® IVC™ Tracer® SX5 Wheelchair

The Invacare IVC Tracer SX5 lightweight frame weighs less than 36 pounds, which makes this chair perfect for rental and long- or short -term use. The Tracer SX5 now combines the design and technology of the 9000 series wheelchairs, making it possible to interchange components. The Tracer SX5 offers multiple seat widths, depths, armstyles, sectional back and seat-to-floor heights by utilizing the dual axle position frame base. The 20" and 22" widths include a semi-pneumatic caster and chrome handrim to accommodate a weight capacity of 300 lb.

#### **FEATURES**

- · Durable, low-maintenance, triple chrome-plated, carbon steel frame is long-lasting
- Dual axle positions allow variation of seat-tofloor heights
- · 14-gauge crossbraces add strength and durability
- Urethane rear tires mounted on "no flex" wheels offer superior performance
- · Sectional back available in 16" and 18" back heights
- · Low-maintenance, dual-sealed precision
- Heavy-duty inner liners keep seat and back from stretching
- · Flame-retardant upholstery resists mildew and bacteria, as well as cold-crack and UV damage
- Meets RESNA standards for manual wheelchairs in product performance

Dimensions	
Seat sizes	
Widths	14", 16", 18", 20", 22"
Depths	16" standard, 18" by
(optional)	seat extension
Seat-to-floor heights	17.5", 19.5" 15.5", 16.5", 18.5" with different wheel and caster packages
Back height	16", 18" sectional
Front riggings	
Hemi swingaway footrest	17"-21"
Hemi swingaway padded elevating legrest	17.5"-21"
Articulating elevating legrest	17"-25"
(Footplates available in aluminum or composite)	
Arm styles Full-length fixed height, padded	
Desk-length fixed height, padded	
Full-length adjustable height, padded	
Desk-length adjustable height, padded	
(All arms are flip back non-removable)	
Weight capacity	250 lb. (300 lb. 20" and 22" width)
Product weight (TRSX5/WD86/28FB)	34 lb.
Limited warranty	lifetime on frame and crossbraces
	5



## Invacare® IVC™ Tracer® SX5 Recliner

The Invacare IVC Tracer SX5 Recliner wheelchair offers the durability of a lightweight manual wheelchair frame with the added versatility and comfort of a full reclining back. The Tracer SX5 Recliner offers dynamic recline range from 100° to 168°. The Tracer SX5 Recliner offers a 300 lb. weight capacity on 20" and 22" widths. This chair also is available with multiple seat widths, depths, armstyles, front riggings, a full range of accessories and a lifetime warranty.

#### **FEATURES**

- The Tracer SX5 Full Recliner offers a dynamic range of recline from 100° to 168°
- Durable, low-maintenance, triple chrome-plated, carbon steel frame
- Urethane rear tires, mounted on "no flex" wheels, offer superior performance
- Low-maintenance, dual-sealed precision bearings
- Heavy-duty inner liners keep seat and back from stretching
- Flame-retardant upholstery resists mildew and bacteria, as well as cold-crack and UV damage
- Upholstery is dual-embossed
- Meets RESNA standards for manual wheelchairs in product performance

Seat sizes	
Widths	14", 16", 18", 20", 22"
Depths	16" standard, 18" by seat extension (optional)
Seat-to-floor height	19.5" 17.5", 18.5" with different wheel and caster packages
Back height	33" recliner (23" high back, 10" removable headrest)
Front riggings	
Hemi swingaway footrest	17"-20"
Hemi swingaway padded elevating legrest	17.5"-21"
Hemi articulating elevating legrest	17"-25"
(Footplates available in aluminum or composite)	
Arm styles Full-length fixed height, padded, removable	
Desk-length fixed height, padded, removable	
Full-length adjustable height, padded, removable	
Desk-length adjustable height, padded, removable	
(All arms are conventional)	
Weight capacity	250 lb. (300 lb. 20" and 22" width)
Product weight (TRSX5/WD86/RC/28)	53 lb.
Limited warranty	lifetime on frame and crossbraces



## Invacare® IVC™ 9000 SL Wheelchair

The Invacare IVC 9000 SL wheelchair is the perfect lightweight rental product for those who want quality and durability at a moderate price. This low-maintenance wheelchair also has a full complement of options and accessories. Its contemporary design and high-performance frame make the 9000 SL a definitively superior value.

#### **FEATURES**

- Durable, low-maintenance, triple chrome-plated, carbon steel frame is long-lasting
- · Multiple frame color variations for customizing
- Dual axle positions allow variations of seat-tofloor heights
- Urethane rear tires mounted on "no flex" wheels offer superior performance
- Urethane casters for better ride and lighter weight
- Low-maintenance, dual-sealed precision
- Heavy-duty inner liners keep seat and back from stretching
- Flame-retardant upholstery resists mildew and bacteria, as well as cold-crack and UV damage
- Space-saver arms reduce overall frame width without limiting seat width
- Meets RESNA standards for manual wheelchairs in product performance

## **Specifications**

Limited warranty

Specifications	
Dimensions	
Seat sizes	
Widths	16", 18", 20"
Depths (optional)	16" standard, 18" by seat extension
Seat-to-floor heights	17.5", 19.5" 15.5",16.5",18.5" with different wheel and caster packages
Back heights - adjustable	15", 16", 17" or 17", 18", 19"
Front riggings	
Hemi swingaway footrest	17"-21"
Hemi swingaway padded elevating legrest	17.5"-21"
(Extra large footplates are also available)	
Articulating elevating legrest	17"-25"
(Footplates available in aluminum or compos	ite)
Arm styles Full-length fixed height, padded, flip back, non-re	movable
Desk-length fixed height, padded, flip back, non-	removable
Full-length adjustable height, padded, flip back, n	on-removable
Desk-length adjustable height, padded, flip back,	non-removable
(All arm styles are space saver)	
Weight capacity	250 lb.
Product weight (9SL/WD86/88)	33 lb.

lifetime on frame and crossbraces



## Invacare® IVC™ 9000 XT Wheelchair

The Invacare IVC 9000 XT wheelchair is a high performance, low maintenance and low total lifetime cost product. The 9000 XT is simply the most versatile lightweight wheelchair available and is the model of choice. This wheelchair boasts a vast offering of frame styles, seat-to-floor heights, back styles, widths, depths, frame colors and a wide range of options and accessories. The 9000 XT components are interchangeable with all other 9000 series and Tracer® series chairs.



9000XT with offset axle and adjustable angle caster will allow additional seat-to-floor heights and increase seating and positioning capabilities. Offset axle offers center-of-gravity adjustment which increases maneuverability and stability.

#### **FEATURES**

- Durable, low-maintenance, chrome-plated, carbon steel frame with aluminum crossbraces
- · Multiple frame color variations for customizing
- Dual axle positions allow variation of seat-to-floor heights
- Offset axle position allows adjustment of center-of-gravity to increase maneuverability or stability to meet each user skill level
- Urethane rear tires mounted on "no flex" wheels offer superior performance
- · Urethane casters for better ride and lighter weight
- Adjustable caster journal can be used with or without offset axle option. (Refer to service manual for a complete list of seat-to-floor heights and angles that are reconfigurable only by the provider.)
- Flame-retardant upholstery resists mildew and bacteria, as well as cold-crack and UV damage
- Meets RESNA standards for manual wheelchairs in product performance



The 9000 XT adjustable-angle back easily changes from -5° to 20° in 5° increments. The back is adjustable from 15" to 19" in 1" increments. Anti-tippers are standard.



The super-hemi (single-axle) seat-to-floor height is achieved with 22" rear wheels and 6" casters. A 13.5" seat-to-floor height is achieved with 20" rear wheels and 6" casters.

Dimensions	
Seat sizes	
Widths	14", 15", 16", 17", 18", 20", 22" adult to hemi frame
	14", 15", 16", 17", 18" super-hemi frame
Depths	16" standard frame
	17", 18" by seat extensions, adult to hemi frame
	16", 17" super-hemi frame (different frame)
Seat-to-floor heights*	17.5", 18.5", 19.5"
	15.5", 16.5" with different wheel and caster packages
<u>-</u>	13.5", 14.5" super-hemi
Back heights	15", 16", 17", 18", 19" adjustable height
	adjustable angle available as option (-5° to 20° adjustment)**
Front riggings	
Hemi swingaway footrest	17"-21"
Hemi swingaway padded elevating legrest	17.5"-21"
Super-hemi swingaway footrest	"- 7"
Super-hemi swingaway padded elevating legrest	"- 7"
Articulating elevating legrest	17"-25"
(Footplates available in aluminum or composite)	
Arm styles	
Full-length fixed height, padded, flip back, removable	
Desk-length fixed height, padded, flip back, removable	
Full-length adjustable height, padded, flip back, remova	ıble
Desk-length adjustable height, padded, flip back, remo	vable
(All arm styles available conventional or space saver)	
Weight capacity	250 lb.
Product weight (9XT/WD86/88)	31 lb.
Limited warranty	lifetime on frame and crossbraces

<sup>\*</sup>See service manual for additional seat-to-floor heights and seat angles available.

<sup>\*\*</sup>Seat angle adjustment of  $0^{\circ}$  and  $6^{\circ}$  that are reconfigurable only by the provider, can be found in the service manual.



#### **FEATURES**

- Durable, low-maintenance, triple chrome-plated, carbon steel frame is long-lasting
- Dynamic recline range from 90° to 168°
- Rear axle base moves, providing stable center of gravity in recline position
- Multiple frame color variations for customizing
- Multiple wheel offerings allow variation of seat-to-floor heights
- Urethane rear tires mounted on "no flex" wheels offer superior performance
- Urethane casters for better ride and lighter weight
- Low-maintenance, dual-sealed precision bearings
- Heavy-duty inner liners keep seat and back from stretching
- Flame-retardant upholstery resists mildew and bacteria, as well as coldcrack and UV damage
- Additional "Lift Packages" assist attendants in positioning patients
- Stocked with attendant (standard) or patient-operated wheel locks
- Meets RESNA standards for manual wheelchairs in product performance

## Invacare® IVC™ 9000 XT Recliner

The Invacare IVC 9000 XT Recliner wheelchair is the number one choice for individuals who demand a truly versatile, yet moderately priced lightweight wheelchair. The 9000 XT Recliner can be configured into an endless number of seat widths, depths and heights, and has multiple frame and upholstery options. The recliner has been standardized to interchange components and accessories with the other 9000 series chairs. Whether you need the lighter recliner package or the assist recliner package, the 9000 XT Recliner can provide it.

Dimensions	
Seat sizes	
Widths	14", 16", 18", 19", 20", 22"
Depths	16", 17", 18"
Seat-to-floor heights	19.5"
· ·	17.5", 18.5" with rear wheel package
Back height	32" recliner (22" high back, 10"
	removable headrest)
Front riggings	
Hemi swingaway footrest	17"-20"
Hemi swingaway padded elevating legrest	17.5"-21"
Hemi articulating elevating legrest	17"-25"
(Footplates available in aluminum or composite)	
Arm styles	
Full-length fixed height, padded, removable	
Desk-length fixed height, padded, removable	
Full-length adjustable height, padded, removable	
Desk-length adjustable height, padded, removable	
(All arms are conventional)	
Weight capacity	250 lb.
Product weight (9RC/WD86/R/28/2R)	51 lb. without lift package
,	62 lb. with lift package
Limited warranty	lifetime on frame and crossbraces



# Invacare® IVC™ 9000Topaz™ Wheelchair

The Invacare IVC 9000 Topaz wheelchair is designed specifically to suit the unique needs of bariatric patients. The Topaz combines the strength and durability of a triple chrome-plated, carbon steel frame with the comfort of heavy-duty front riggings and arms for proper positioning. With a variety of seat widths and depths, frame reinforcements that include gussetted sides, caster journals and dual heavy-duty crossbraces, this chair provides not only ease in maneuverability but stability that will be trusted for years to come.

#### **FEATURES**

- · Available in vinyl or nylon upholstery
- Durable low-maintenance, triple chrome-plated, carbon steel frame is long-lasting
- Reinforced frame includes gussetted sides, caster journals and dual heavy-duty crossbraces
- Multiple frame color variations available for customizing
- Heavy-duty front riggings provide proper positioning for feet with footrest and elevating legrest
- Heavy-duty arms designed to accommodate weight capacity and provide proper positioning
- Dual axle positions allow for quick and simple seat-to-floor height variations, 17.5" and 19.5"
- Low-maintenance, dual-sealed precision bearings
- Heavy-duty triple inner liners keep seat and back from stretching
- Flame-retardant upholstery resists mildew and bacteria, as well as cold-crack and UV damage

Dimensions	
Seat sizes	
Widths	20", 22", 24", 26", 28", 30'
Depths	18", 20"
Seat-to-floor heights	17.5", 19.5" 15.5" with different wheel and caster packages
Back heights - adjustable	16", 17", 18"
Front riggings	
Swingaway footrest	
Swingaway padded elevating legrest	
Footrest with one piece wooden footboard	
Arm styles Full-length fixed height, padded, removable	
Desk-length fixed height, padded, removable	
Full-length adjustable height, padded, removable	
Desk-length adjustable height, padded, removable	
(All arm styles available conventional)	
Weight capacity	700 lb.
Product weight (9TPZ/WD2818/28)	85 lb.
Limited warranty	lifetime on frame and crossbraces



# Invacare® IVC™ 9000 XDT Wheelchair

The Invacare IVC 9000 XDT wheelchair is designed for individuals who demand an extra wide, heavy-duty or tall chair with seat-to-floor heights ranging from hemi through tall. The high-performance 9000 XDT provides a 350-lb. weight capacity and components interchangeable with the other 9000 series wheelchairs.

#### **FEATURES**

- Durable, low-maintenance, triple chromeplated, carbon steel frame is long-lasting
- Tall frames are designed to accommodate taller individuals and provide a more comfortable ride
- · Multiple frame color variations for customizing
- Triple-axle positions allow wide variation of seat-to-floor heights
- Urethane rear tires mounted on "no flex" wheels offer superior performance
- Urethane casters for better ride and lighter weight
- Low-maintenance, dual-sealed precision bearings
- Heavy-duty inner liners keep seat and back from stretching
- Flame-retardant upholstery resists mildew and bacteria, as well as cold-crack and UV damage

Dimensions	
Seat sizes	
Widths	16", 17", 18", 19", 20", 22" tall to adult to
	hemi frame
Depths	16", 17", 18", 19", 20"
Seat-to-floor heights	17.5", 19.5" 21.5" option
Back heights - adjustable	15", 16", 17", 18", 19"
Front riggings	
Hemi swingaway footrest	17"-21"
Hemi swingaway padded elevating legrest	17.5"-21"
Super-hemi swingaway footrest	11"-17"
Super-hemi swingaway padded elevating legrest	11"-17"
Articulating elevating legrest	17"-25"
(Footplates available in aluminum or composite)	
Arm styles Full-length fixed height, padded, flip back, removable	
Desk-length fixed height, padded, flip back, removable	е
Full-length adjustable height, padded, flip back, remov	rable
Desk-length adjustable height, padded, flip back, remo	ovable
(All arm styles available conventional or space saver)	
Weight capacity	350 lb.
Product weight (9XDT/WD08/28)	37 lb.
Limited warranty	lifetime on frame and crossbraces





Heavy-duty package increases weight capacity to 450 lb. and consists of heavy-duty composite wheels with pneumatic tires with flat-free inserts and semi-pneumatic casters.

## Invacare® IVC™ Tracer® IV Wheelchair

The Invacare IVC Tracer IV wheelchair features superior durability, rollability and streamlined looks. It's designed for the active individual who needs the strength of a reinforced heavyduty frame. Its many options include components interchangeable within the Tracer series. The Tracer IV is long-lasting and low-maintenance.

#### **FEATURES**

- Dual axle position allows variation of seat-tofloor heights
- Durable, low-maintenance, triple chromeplated, carbon steel frame is long-lasting
- Reinforced frame includes gussetted sides, caster journals and dual crossbraces
- Corded urethane rear tires mounted on "no flex" wheels offer superior performance (350 lb.)
- Low-maintenance, dual-sealed precision bearings
- Heavy-duty inner liners keep seat and back from stretching
- Flame-retardant upholstery resists mildew and bacteria, as well as cold-crack and UV damage
- · Upholstery is dual-embossed

Dimensions	
Seat sizes	
Widths	18", 20", 22", 24"
Depth	18" (standard),
0	20" (by seat extension)
Seat-to-floor height	17.5", 19.5"
Back height	16" fixed
Front riggings	
Swingaway footrest	17"-21"
Swingaway padded elevating legrest	17.5"-21"
Articulating elevating legrest	17"-25"
(Footplates available in aluminum and composite)	
Arm styles	
Full-length fixed height, padded, removable	
Desk-length fixed height, padded, removable	
Full-length adjustable height, padded, removable	
Desk-length adjustable height, padded, removable	
(All arms are conventional)	
Weight capacity	350 lb. and 450 lb.
Product weight (T4/WD28/28)	42 lb.
Limited warranty	lifetime on frame and crossbraces



## Invacare® IVC™ 900 Wheelchair

The Invacare IVC 900 series wheelchair is designed to meet the physical needs of the most demanding consumer. A reinforced 900 chair leads the way when it comes to durability and reliability: we have strengthened the 900 in numerous ways, including additional supports, gussets and reinforcements. The 900 redefines the term "custom wheelchair". Each one is built to the exact specification of the consumer. The chair is available in multiple widths and depths, back styles, seat-to-floor heights, and weight limits up to 350 lb., resulting in a chair that can accommodate almost any user.

## **FEATURES**

- · Steel frame with durable finish
- Multiple width and depth combinations can accommodate almost any consumer
- · Gusseted frame option withstand heavy use
- Dual crossbrace option provides extra strength and rigidity
- Shock absorbing forks handle the abuse the front-end may receive

4", 15", 16", 17", 18", 3", 20", 21", 22", 23", 24" 3", 14", 15", 16", 17",
9", 20", 21", 22", 23", 24" 3", 14", 15", 16", 17",
3", 19", 20", 21", 22"
3" hemi, 20" adult 9", 21", 22", 23", 24" v frame modification
4"-24" sectional back 4" -24" permanent ked back height " increments)
50 lb., 275 lb., 00 lb., 350 lb.
ependent on custom onfigurations
years on frame and ossbraces



# Invacare® IVC™ 900 Recliner

The Invacare IVC 900 bariatric recliner is custom-made with multiple widths and depths while offering a 350 lb. weight capacity. The 900 recliner, like the IVC 900 wheelchair, leads the way in durability and reliability. Frame gussets and reinforcements allow this chair to meet the special needs of those individuals who need a reclining back wheelchair that accommodates a weight capacity of 350 lb.

## **FEATURES**

- · Steel frame with durable finish
- Multiple width and depth combinations can accommodate almost any user
- · Gusseted frame option withstands heavy use
- Dual crossbrace option provides extra strength and rigidity
- Shock absorbing forks handle the abuse the front-end may receive

Dimensions					
Seat sizes					
Widths	14", 15", 16", 17", 18", 19", 20", 21", 22", 23", 24"				
Depth	13", 14", 15", 16", 17", 18", 19", 20", 21", 22"				
Seat-to-floor height	18" hemi, 20" adult 19", 21", 22", 23", 24" by frame modification				
Back height	32" reclining 22" high back and 10" removable headrest				
Arm styles					
Full-length fixed height, padded, removable					
Desk-length fixed height, padded, removable					
Full-length adjustable height, padded, removable					
Desk-length adjustable height, padded, removable					
(All arms are conventional)					
Weight capacity	350 lb.				
Product weight	dependent on custom configuration				
Limited warranty	5 years on frame and crossbraces				



# Invacare® IVC™ 9000 Jymni™ Wheelchair

The Invacare IVC 9000 Jymni wheelchair offers your growing list of pediatric clients versatility, functionality and maneuverability. These chairs are specifically designed to suit the needs of your smaller clients and promote many outstanding features.

#### **FEATURES**

- Durable, low-maintenance, chrome-plated, carbon steel frame is long-lasting
- Stroller handles provide ergonomic positioning for caregiver and remove for ease in transporting
- Adjustable front riggings provide ideal positioning for pediatric clients
- Urethane rear tires mounted on "no flex" wheels offer superior performance
- Virtually no-maintenance, dual-sealed precision bearings
- Flame-retardant upholstery resists mildew and bacteria, as well as cold-crack and UV damage
- Dual-axle positions allow for quick and simple seat-to-floor height adjustments of 15" and 16"
- Seat width and depth options provide proper positioning
- Adjustable angle back easily changes from -5° to 20° in 5° increments
- Meets RESNA standards for manual wheelchairs in product performance

Dimensions	
Seat sizes	
Widths	10", 12", 14", 16"
Depths	10", 12", 14"
Seat-to-floor heights	15", 16" 17" with different wheel and caster packages
Back heights - adjustable	14"-15"
Front riggings	
70° swingaway footrest	7"-11"
80° swingaway footrest	11"-15"
90° swingaway footrest	5"-10"
Elevating legrest	7"-14" or 11"-15"
Arm styles Full-length fixed height, padded, flip back, re	emovable
Desk-length fixed height, padded, flip back,	removable
Full-length adjustable height, padded, flip ba	ck, removable
Desk-length adjustable height, padded, flip	back, removable
(All arm styles available conventional or sp	ace saver)
Weight capacity	125 lb.
Product weight (9JYLT/WD22/88JY)	28 lb.
Limited warranty	5 years on frame and crossbraces





9000 Jymni Recliner allows for a dynamic recline range of 95° to 160°. Multiple frame color variations for customizing.

# Invacare® IVC™ 9000 Jymni™ Recliner Wheelchair

The Invacare IVC 9000 Jymni Recliner wheelchair offers your growing list of pediatric clients versatility, functionality and maneuverability. These chairs are specifically designed to suit the needs of your smaller clients and provide many outstanding features.

#### **FEATURES**

- Durable, low-maintenance, chrome-plated, carbon steel frame is long-lasting
- Dynamic recline from 95° to 160°
- Stroller handles provide ergonomic positioning for caregiver and remove for ease in transporting
- Adjustable front riggings provide ideal positioning for pediatric clients
- Flame-retardant upholstery resists mildew and bacteria, as well as cold-crack and UV damage
- Seat width and depth options provide proper positioning
- Meets RESNA standards for manual wheelchairs in product performance

Dimensions Seat sizes	
Widths	10", 12", 14", 16"
Depths	10", 12", 14"
Seat-to-floor heights	15", 16" 17" with different wheel and caster packages
Back heights	27" reclining 16" high back and 11" removable headrest
Front riggings	
70° swingaway footrest	7"-11"
80° swingaway footrest	11"-15"
90° swingaway footrest	5"-10"
Elevating legrest	7"-14" or 11"-15"
Arm styles Full-length fixed height, padded, flip back, r	emovable
Desk-length fixed height, padded, flip back	removable
Full-length adjustable height, padded, flip b	ack, removable
Desk-length adjustable height, padded, flip	back, removable
(All arm styles available conventional)	
Veight capacity	125 lb.
roduct weight (9JYLT/WD42/28JY/R)	41 lb.
imited warranty	5 years on frame and crossbraces



Model no. ICR330

# Invacare® Ramps\*

Designed for ease of use and durability, these folding ramps offer an economical alternative to the motorized lift. They're ideal for both manual and power wheelchairs, as well as 3- or 4-wheel scooters. Plus, they're extremely lightweight for effortless transport. Clearly, these ramps are a cost-effective way to add to your Total One Stop Shopping⁵ purchases with Invacare.

#### **FEATURES**

- · Sturdy and dependable 600-lb. weight capacity
- Ideal for manual or power wheelchairs and 3or 4-wheel scooters
- Convenient carrying handle is standard on all ramps
- Heavy-duty industrial Velcro® strap ensures ramp security during transportation
- Transports and stores like a suitcase
- High-traction tape on the surface is long-lasting and maintenance-free
- · Stocked in Invacare's warehouse



Unique folding design makes ramps easy to carry.

## **Specifications**

#### Dimensions

	Model no. ICR630	Model no. ICR330
Size	6' l x 30" w	3' l x 30" w
Style	multifold	single fold
Folded size	40.625" l x 16" w x 8.5" d	36" l x 16" w x 4" d
Product weight	28 lb.	14 lb.
Weight capacity	600 lb.	600 lb.
Limited warranty	l year	l year

<sup>\*</sup>Ramps are manufactured by Prarie View Industries and distributed by Invacare Corporation.



# Invacare® IVC™ CLD

The Invacare IVC CLD (Cyclical Lever Drive) was designed for individuals to self-propel themselves with one arm. The opportunity to be more active and mobile which allows individuals to enjoy their independence - and the use of the CLD supports this. With its front-caster steering mechanism, simple rowing motion design and adjustability in height and stroke length, Invacare's CLD is the answer to an easy-to-operate one-arm drive.



Non-removable Device (FR)

1820 9000 XT, XT Recliner, SL, XDT, Jymni, Jymni Recliner, Tracer SX5, SX5 Recliner, Tracer IV, EX2, 900, 900 Recliner


1335 (14"-22" w) - Tracer EX2, IV,



Non-folding Device

Axles

# I. 1577QR - Quick Release -9000 Jymni, 9000 XT, SL, Tracer



- 3.  $9357L-Amputee\ Attachment$ 9000 XDT, XT, SL, Tracer SX5, EX2
- 4. OFFAXL Off axle 9000 XT

## **Specifications**

9.1 lb.
1.5"
15"-16"
15"-18"
black/chrome
right- or left-hand
250 lb.
Two years*

\*Invacare warrants all Cyclical Lever Drive mechanisms to be free from defects of materials and workmanship for a period of two years from date of purchase. See warranty for additional details.

There are no service requirements for the CLD, although service requirements for the

Kits available for provider installation. Right- or left-handed configuration. Provider installed kits for Tracer EX2, Tracer SX5 and 9000 SL. See 9000 XT price list (95-66 Rev. 10/02) for availability.



Wheel lock position.



Easy to disengage.



Wheel braking position.



Overhead Anti-Theft Device

- 1514 Tracer EX2, SX5, Invacare Transport, SX5 Recliner, Tracer IV, (18"&20")
- 1715 Single-Pole (not shown) Tracer IV



Accessories

- I. 1300 Utility Tray 9000 XDT, XT, XT Recliner, 900, 900 Recliner
- 2. I209 Tire Pump XT Recliner, 900, 900 Recliner
- I305L Crutch & Cane Carrier – 9000 Topaz<sup>™</sup>, XDT, XT.Tracer IV
- 4. 1100 Tool Kit 9000 XDT, XT, SL, XT Recliner
- 5. I571 Fabric Clothing Protectors – 9000 XDT, XT




Upholstery

- ImpactG Impact Guards 9000 XT
- STC Stencil Back Upholstery 9000 XDT, XT, XT Recliner, Tracer IV, SX5, EX2, SX5 Recliner
- 3. SP-01 Seat Pouch 9000 XT
- 4. BP-02 Backpack 9000 XT (shown in Large 16" and Small 14"-16")



Seat/Chest Positioning Straps

- I. 1321BK/7321BK Hook and Loop - not available on 9000 Topaz,Transport
- 2. 1515 Auto Style Buckle not available on Transport
- 3. 1311BK/7311BK Airline Buckle - 9000 Jymni, XT, SL, XDT



Upholstery

- 9130 Solid Back Insert 9000 XDT, XT, Jymni, Jymni Recliner, 900 (T9130 - 9000 Topaz)
- 2. 9662L Solid Seat Insert Tracer SX5, SX5 Recliner, EX2, 9000 XT, XT Recliner, XDT
- 3. 9120 Solid Seat Insert 9000 SL, XT, XDT, XT Recliner, 9000 Jymni, Jymni Recliner, Tracer IV, 900, 900 Recliner (T9120 - Topaz)
- 1608 Solid Back Insert 9000 XT Recliner, 900 Recliner



Anti-Tippers

- 9758 Rear Anti-Tippers EX2 (FF only)
- 2. 1360SH Super Low Adjustable – 9000 XT
- 3. I360A Adjustable Rear Anti-Tippers – 9000 XDT, XT, SL. Tracer EX2, SX5
- 4. 1360 Adjustable Rear Anti-Tippers – 9000 XT, SL, Jymni, Tracer EX2, SX5, Tracer IV
- 5. 2226 Front Anti-Tippers 9000 XDT, XT, SL, XT Recliner



Tires

- I. Urethane
- 2. Urethane High Profile
- 3. Pneumatic Flat-Free All tires available in 20", 22" and 24" sizes

Available on all chairs



IV/O<sub>2</sub> Holders

- 1. 1497 O<sub>2</sub> Holder Package Assembly – 9000 XDT, XT, SL, Tracer IV, SX5, EX2 (T1497 -Topaz, 1488 - 900, 900 Recliner, 9000 Jymni, Jymni Recliner)
- 2. 1496 IV/O<sub>2</sub> Holder Package Assembly – 9000 XDT, XT, SL, Tracer IV, SX5, EX2 (T 1496 -9000 Topaz) (1487 - 900, 900 Recliner, 9000 Jymni, Jymni Recliner)
- 3. 2218 IV Rod, Telescopic 9000 XT, SL, XDT, Tracer IV, EX2, SX5
- 4. 2219 IV Rod, Telescopic XT Recliner
   1322SXR – IV Rod and Bracket – Tracer SX5 Recliner (not shown)



Handrims

- I. CHDRM Chrome-Plated 9000 XDT, Topaz, XT, SL, Jymni Tracer IV, SX5, EX2, SX5 Recliner
- 2. COM Composite 9000 XDT, XT, SL, Tracer SX5, EX2, SX5 Recliner
- 3. AHDRM Aluminum 9000 XDT, XT, SL, SX5, SX5 Recliner
- 4. 1512 Plastic Coated 9000 XDT, 9000 XT
- 5. 1225-8 8 Aluminum Vertical Projection Handrim Plastic Coated – 9000 XDT, 9000 XT
- 1225-12 12 Aluminum Vertical Projection Handrim Plastic Coated – 9000 XDT, XT, XT Recliner
- 1224-8 8 Aluminum Oblique Projection Handrim Plastic Coated – 9000 XT Recliner
   BHDRM – Black Handrim - 9000 XDT, Jymni (not shown)



Casters

- 1. 1292 8"x 1.25" Pneumatic 9000 XT, SL
- 2. AT903 8"×1.25" Composite Urethane – 9000 XDT, XT, SL, Tracer® SX5, SX5 Recliner, Tracer IV
- 3. 1255 8"x1" Solid Rubber Tracer SX5, EX2
- 4. AT902 6" Composite Urethane, 3 Spoke – 9000 XT, SL, Jymni™, Tracer EX2, SX5
- 1252 8"x 1.25" Pneumatic Tracer SX5, EX2, SX5 Recliner
- 6. 1228 8"x 2" Semi-Pneumatic 9000 XDT, Topaz<sup>™</sup>, 9000 XT, Tracer IV, SX5, SX5 Recliner
- 7. 1687 6"x1" Dark grey composite Urethane 9000 Jymni
- 8. 1259 6"x2" Semi-Pneumatic 9000 Topaz
- 9. 1290 8"x I" Solid Rubber, Charcoal – 9000 XDT, XT, SL, Tracer IV
- 10.4252 8"x1.25" Pneumatic with flat-fee inserts - Tracer EX2, SX5, SX5 Recliner
- $$\label{eq:local_local_state} \begin{split} &11.\,1492-8"x\,1.25"\,\,Pneumatic\\ &\text{with flat-fee inserts}-\,\,9000\,\,XT,\\ &9000\,\,SL \end{split}$$





Wheel Locks



Conventional Arms

- 9690 Pull to Lock –9000 XDT, Topaz, 9000 XT, SL, Tracer IV 9691 – Push to Lock (Hemi) – 9000 XDT,Topaz, XT, SL, Jymni, Tracer IV
- I. 9697R Pull to Lock Hill Holder – 9000 XT Recliner
- 2. 9695 Hill Holder 9000 XT, SL, 900, 900 Recliner
- 3. A1567 Pull to Lock Clamp-on 9000 XT
- 4. 1721 High Push to Lock Clamp-on - 9000 XT
- 5. I 329BK Toggle Wheel Lock Extension – 9000 XT Recliner, Tracer IV, EX2, SX5, SX5 Recliner, 900, 900 Recliner
- 18 Full-Length Fixed Height Tracer IV, EX2, Transport, SX5 Recliner, 9000 XT, XDT, XT Recliner (18JY - Jymni, 18FB - Flip-back - SX5)
- 19 Full-Length Adjustable Height – Tracer IV, EX2, SX5 Recliner , 9000 XT, XDT, XT Recliner (19]Y - Jymni, XT Recliner, 19FB - Flip-back - SX5)
- 3. 28 Desk-Length Fixed Height – Tracer IV, EX2, SX5 Recliner, 9000 XT, XDT, XT Recliner (28JY - Jymni, 28FB - Flip-back -SX5)
- 4. 29 Desk-Length Adjustable Height – Tracer IV, EX2, SX5 Recliner, 9000 XT, XDT, XT Recliner (29JY - Jymni, XT Recliner, 18FB - Flip-back - SX5)
- TZ33A Swingaway footrest 9000 Topaz
- 2. TZ34A Swingaway elevating legrest 9000 Topaz
- TZ1263 Footrest with onepiece wooden footboard -9000 Topaz



9000 Topaz Arms

- I. T18 Full-Length Fixed Height
- T19 Full-Length Adjustable
  Height
- 3. T28 Desk-Length Fixed Height
- 4. T29 Desk-Length Adjustable Height

Available on 9000 Topaz.



9000 Topaz Front Riggings



Lightweight Space Saver Arms

- 1. 79 Full-Length Adjustable Height – 9000 XDT, XT, SL, 79JY - 9000 XT, Jymni
- 2. 78 Full-Length Fixed Height 9000 XDT, XT, SL 78JY 9000 Jymni
- 3. 89 Desk-Length Adjustable Height – 9000 XDT, XT, SL 89JY - 9000 XT, Jymni
- 4. 88 Desk-Length Fixed Height – 9000 XDT, XT, SL 88|Y - 9000 Jymni



9000 Jymni Front Riggings

- I.  $J932A 70^{\circ}$  swingaway footrest
- 2. J932A-14 80° swingaway footrest
- 3.  $J933A 90^{\circ}$  swingaway footrest
- 4. J934A elevating legrest
- 5. J934A-14 elevating legrest

Available on 9000 Jymni

All have aluminum footplates



Footrests

- T963A/C 60° Elevating
   Footrest Aluminum/Composite
- 2. T93HA/C/S Aluminum/ Composite -Tracer EX2, SX5, SX5 Recliner, 9000 SL, XT, XT Recliner, XDT, Tracer IV
- 3. T93HE Composite Tracer EX2, SX5



**Footrests** 

- SH93F Superlow Extension Swingaway, Aluminum – 9000 XT
- 2. H93M Hemi Minimum Swingaway, Aluminum – 9000 XT, XDT, XT Recliner
- 3. H93E Hemi Extension Swingaway, Aluminum – 9000 XT, XDT, XT Recliner
- SH93ES Super Low Extension Swingaway, Aluminum – 9000 XT



Legrests

- I. AHS4A Hemi Smartleg Articulating Legrest, Aluminum – 9000 XT, XT Recliner, Tracer SX5, SX5 Recliner, Tracer IV SH904M – Superlow Elevating Legrest, Aluminum – 9000 XT
- T94HA/C Hemi Elevating Legrest Aluminum/Composite – Tracer SX5, SX5 Recliner, EX2, 9000 SL, XT, XT Recliner XDT, Tracer IV

(not shown)

- T94HS Hemi Elevating Legrest Aluminum 14" – Tracer SX5, SX5 Recliner, 9000 XT, XT Recliner, Tracer IV (not shown)
- 3. SH904MS Superlow Elevating Legrest, Aluminum – 9000 XT
- 4. AHL4A/C Hemi Articulating Legrest, Aluminum/Composite – 9000 XDT, XT, SL, XT Recliner, Tracer EX2, SX5, SX5 Recliner, Tracer IV



Straps

- I. 1346 Heel Loops 9000 XDT, XT, SL, XT Recliner, Tracer IV, SX5, EX2, SX5 Recliner
- 2. 1345 Toe Loops XT Recliner
- 3. 1600 BK Heel Loops with Ankle Straps – 9000 XDT, XT, SL, XT Recliner, Tracer IV
- 4. 1348 Legrest Panel XT Recliner
- 1337 Calf Strap 9000 XDT, XT, SL, XT Recliner
- 6. 9342 H Strap 9000 XDT, XT, XT Recliner, Tracer IV



Footplates

- 1. 1350 Extra Large Aluminum 9000 Topaz™ (standard), 9000 XDT, XT, SL, XT Recliner, Tracer SX5, EX2, IV, SX5 Recliner, Tracer IV
- AT5543 Adjustable Angle 9000 XDT, XT, XT Recliner, Tracer SX5 Recliner, EX2
- C Composite 9000 XDT, XT, SL, XT Recliner, Tracer SX5, EX2, SX5 Recliner, Tracer IV
- 4. A Aluminum 9000 XDT, XT, SL, XT Recliner, Tracer SX5, EX2, SX5 Recliner, Tracer IV

N 4		. 1	
I۲I	വ	eı	

		Tracer EX2	Transport	Tracer SX5	Tracer SX5 Recliner	9000 SL	9000 XT	9000 XT Recliner	9000 Topaz	9000 XDT	Tracer IV	900	900 Recliner	9000 Jymni	9000 Jymni Recliner
Description	Part no.												,		
Non-removable Arms w/ Screws	NRAS														
Non-removable Device (FR)	1820														
One-Arm Drive				*											
One-Piece Wood Footboard	T1263														
O <sub>2</sub> Holder	1497								T1497						
O <sub>2</sub> Holder	1488														
Overhead Anti-Theft Device	1514										18", 20" only				
Patient Operated Wheel Locks (Recliner)	1245														
Pull-to-Lock Clamp-On Wheel Locks	A1567														
Pull-to-Lock Wheel Lock	9690														
Pull-to-Lock Hill Holder	9697R														
Push-to-lock Wheel Lock	9691														
Push-to-lock Wheel Lock (Red	c.) 9696R														
Push Bar	T1483														
Quick Release Axles	1577QR														
Rear Anti-Tippers	1358M														
Rear Anti-Tippers	T1358														
Rear Anti-Tippers	1360														
Rear Anti-Tippers	1360A														
Rear Anti-Tippers	1360B														
Rear Anti-Tippers	9758 I	FF only													
Seat Pouch	SP-01														
Solid Back Insert	9130								T9130						
Solid Back Insert (Recliner)	1608														
Solid Seat Insert	9120								T9120						
Stencil Back Upholstery	STC														
Stroller Handle	1483JY														
Super-Low Adj. Rear Anti-Tipper	1360SH														
Super-Low Amputee Attach.	9357SH														
Telescopic IV Rod	1322SXR														
Telescopic IV Rod and Bracke	tT1322														
Telescopic IV Rod	2218														
Telescopic IV Rod (Recliner)	2219														
Tire Pump	1209														
Toggle Tool Kit	1100														
Torso Support (Recliner)	3040											1164			
Utility Tray	1300														
Wheel Lock Extensions	1329BK														
Wheel Lock Extensions	A1576														
*CLD kit.															

\*CLD kit.

								Mo	del						
		Tracer EX2	Transport	Tracer SX5	Tracer SX5 Recliner	9000 SL	9000 XT	9000 XT Recliner	9000 Topaz	9000 XDT	Tracer IV	900	900 Recliner	9000 Jymni	9000 Jymni Recliner
Description	Part no.	-						-							
Adjustable Angle Footplate	AT5543														
Adjustable Caster Journal	ADJC														
Adjustable Tension Back Upholstery	ADJTENS														
Amputee Attachment	9357L											9357	9357		
Auto-style Buckle Seat Positioning Strap	1515	•							T1515						
Auto-style Buckle Chest Positioning Strap	9515						•	-					-	•	-
Airline Buckle Seat Positioning Strap	1311BK														
Airline Buckle															
Chest Positioning Strap	7311BK														
Backpack	BP02														
Calf Strap	1337														
Carry Pocket	1160														
Carrying Straps	1409														
Crutch & Cane Carrier	1305L											1305	1305		
Drop Hook Solid Folding Seat	9662L														
Dual Back Liner	RB														
Dual Seat Liner	RS														
Extra Large Footplates	1350								std						
Fabric Clothing Protectors	1571														
Footplate Parallel to floor	FPP														
I" Forward Positioner	9357JY														
Front Anti-Tippers	2226											1326	1326		
"H" Strap	9342														
Headrest w/Push Handle ext.	1483R														
Head Support (Recliner)	3046											3045			
Heel Loops	1346														
Heel Loops w/Ankle Strap	1600BK														
High Push-to-lock clamp-on wheel locks	1721														
Hill Holder	9695														
Hook and Loop Seat Positioning Strap	1321BK														
Hook and Loop Chest Positioning Strap	7321BK												-	7321	7321
Impact Guards	IMPACTG														
IV/O <sub>2</sub> Holder	1487														
IV/O <sub>2</sub> Holder	1496								T1496						
Legrest Panel	1348														
Minimum Legrest Ext.	MFXL														
Nonfolding Device	1335														













	(A) (B) (B)		0.0		C (8 -772			A 6-20		A 6 472		0 0 0	
	Tracer® EX	2 MANU	Transport		racer® SX!	;	Tr	acer® SX! Recliner	5	9000 SL	центу	9000 XT VEIGHT	
HCPC Codes assigned Eligible for no charge upgrade <sup>†</sup>	K0002	MAING	E1038		.0003, K000 (0001,K000 K0028 <sup>7</sup>			003 & K00 0006 &K00		K0004		K0004 0001,K0002	2,K0003
Widths	16", 18", 20	"	18", 20"		14", 16", 18" 20", 22"	,	ı	4", 16", 18" 20", 22"	',	16", 18", 20	)" [4	4" - 18", 20' 1" increme	*
Depths	16" 18" <sup>2</sup>		16"		16" 18" <sup>2</sup>			16" 18" <sup>2</sup>		16" 18" <sup>2</sup>		16" 17"², 18"	2
FF	23.5", 25.5", :: 22.75", 24. ble: 24.5", 26	75"	20.5", 22.5'	' 22	.5", 24.5", 26 28.5", 30.5"	5.5"		5", 24.5", 26 28.5", 30.5"		3.25" <sup>6</sup> , 25.2 27.25" <sup>6</sup>		25", 22.25", 24.25",  25.2 27.25", 29.2	25",
Overall depth	32"		22"		32"		36"	(53.5" in rec	line)	30"		30"	
Seat-to-floor heights	17.5" 19.5" 5",³ 16.5",³18		19.5"	15.!	17.5" 19.5" 5",³16.5",³ 18	3.5"³	ı	19.5" 7.5",³ 18.5"	' <sup>3</sup> 15	17.5" 19.5 5.5",³ 16.5",³	18.5"3	7.5", 18.5", 15.5", 16. er-low: 13.5	5'
Back styles "increments	fixed 16"		9" fold dow 9.5" (folded		16", 18" sectional		(2	eclining 33 3" high bac headrest e	:k,	adjustable 15"-17" 17"-19"		adjustable	
Chair weight w/o riggings	36 lb.		27 lb.		34 lb.			53 lb.		33 lb.		31 lb.	
Weight capacity	250 lb.		250 lb.		250 lb. 300 lb. <sup>5</sup>			250 lb. 300 lb. <sup>5</sup>		250 lb.		250 lb.	
Warranty frame/crossbraces	lifetime		3 years		lifetime			lifetime		lifetime		lifetime	
Arm type 08 - Permanent	offset padded		padded		n/a			n/a		n/a		n/a	
18 - Full-length fixed height, conventional	padded removable	!	n/a		dded, flip ba on-removab 18FB			padded removable		n/a	Р	added, flip b removabl	
28 - Desk-length fixed height, conventional	padded removable	· ·	n/a		dded, flip ba on-removab 28FB			padded removable		n/a	P	added, flip b removabl	
19 - Full-length adjustable height, conventional	padded removable		n/a		dded, flip ba on-removab 19FB			padded removable		n/a	P	added, flip b removabl	
29 - Desk-length adjustable height, conventional	e padded removable	!	n/a		dded, flip ba on-removab 29FB			padded removable		n/a	Р	added, flip t removabl	
78 - Full-length fixed height, space saver	n/a		n/a		n/a			n/a		dded, flip b on-remova		added, flip t removabl	
88 - Desk-length fixed height, space saver	n/a		n/a		n/a			n/a		dded, flip b on-remova		added, flip l removabl	
79 - Full-length adjustable height, space saver 79]Y also available on 9000 XT	n/a		n/a		n/a			n/a		dded, flip b		added, flip t removabl	
89 - Desk-length adjustable height, space saver 89JY also available on 9000 XT	e n/a		n/a		n/a			n/a		dded, flip b on-remova		added, flip b removabl	

<sup>†</sup> See Program Memorandum (Transmittal B-01-68) for policies and instructions. 1. 20" and 22" only. 2. By seat extension. 3. Available with different wheel and caster packages. 4. By frame modification

















	. O Common of the common of th	O.A.E.	A 8 507	- 4 E		22	<b>A</b>	Of the state of th
	9000 XT Recliner	9000 Topaz™	9000 XDT	Tracer® IV	900	900 Full Recliner	9000 Jymni™	9000 Jymni™ Recliner EDIATRIC
K	0004 & K0028	K0007	K0007	K0007	K0001 & K0006	K0001 & K0028 K0006 & K0028	K0009	K0009 & K0028
	14", 16", 18", 19", 20", 22"	20", 22", 24" 26", 28", 30"	16-20", 22" I" increments	18", 20", 22", 24"	I4"-24" I" increments	I4"-24" I" increments	10", 12", 14", 16"	10", 12", 14", 16"
	16" 17"², 18"²	18",20"	16"-20" I" increments	18", 20" <sup>2</sup>	13"-22" I" increments	13"-22" I" increments	10", 12", 14"	10", 12", 14"
	.5", 23.5", 25.5" .5", 27.5", 29.5"	27.5", 29.5", 31.5" 33.5", 35.5", 37.5"	23.5", 24.5" 25.5" 26.5", 27.5", 29.5"	26.25", 28.25", 30.25", 32.25"	21.5"-31.5" I" increments	21.5"-31.5" 1" increments	17.5", 19.5", 21.5"	17.5", 19.5", 21.5"
36"	(53.5" in recline)	35"	30"	33.5"	32"	32"	27.5"	27.5"
	19.5" 17.5",³ 18.5"³	17.5", 19.5", 15.5"³	17.5", 19.5", 21.5" <sup>4</sup>	17.5", 19.5"	18" hemi 20" adult (19", 21", 22", 23", 24"	18" hemi 20" adult )4 (19", 21", 22", 23",	15", 16", 17" <sup>3</sup> 24") <sup>4</sup>	15", 16", 17" <sup>3</sup>
(2	reclining 33" 23" high back, ' headrest ext)	adjustable 16", 17", 18"	adjustable 15"-19"	fixed 16"	sectional 14"-24" permanent 14"-24"	reclining 32" (22" high back, 10" headrest ext.)	adjustable 14"-15"	reclining 27" (16" high back, 11" headrest ext.)
	lb. w/o lift pkg. 2 lb. w/lift pkg.	85 lb.	37 lb.	42 lb.	differs per configuration	differs per configuration	28 lb.	41 lb.
	250 lb.	700 lb.	350 lb.	350 lb. 450 lb. <sup>3</sup>	250 lb./275 lb. 300 lb./350 lb.	350 lb.	125 lb.	125 lb.
	lifetime	lifetime	lifetime	lifetime	5 years	5 years	5 years	5 years
	n/a	n/a	n/a	n/a	n/a	n/a	n/a	n/a
	padded removable	padded removable T18	padded, flip back removable	padded removable	padded removable	padded removable	padded, flip back removable 18JY	padded removable 18JY
	padded removable	padded removable T28	padded, flip back removable	padded removable	padded removable	padded removable	padded, flip back removable <sup>28J</sup> Y	padded removable <sup>28</sup> JY
	padded removable	padded removable	padded, flip back removable	padded removable	padded removable	padded removable	padded, flip back removable 19JY	padded removable 19JY
	padded removable	padded removable T29	padded, flip back removable	padded removable	padded removable	padded removable	padded, flip back removable 29 Y	padded removable 28JY
	n/a	n/a	padded, flip back removable	n/a	padded removable	n/a	padded, flip back removable 79JY	n/a
	n/a	n/a	padded, flip back removable	n/a	padded removable	n/a	padded, flip back removable 88JY	n/a
	n/a	n/a	padded, flip back removable	n/a	padded removable	n/a	padded, flip back removable 79JY	n/a
	n/a	n/a	padded, flip back removable	n/a	padded removable	n/a	padded, flip back removable 89JY	n/a

dification. 5. 300 lb weight capacity on 20" and 22" widths. 6. Add 3/s" to overall width for space saver adjustable arms, add 1/4" to overall width for conventional arms. 7. Recliner retro kit available.

Model Number	Description
Transport	
TRAN86	16" seat width, Permanent Full Padded Arm, Fold Down Back, Black Nylon
TRAN06	18" seat width, Permanent Full Padded Arm, Fold Down Back, Black Nylon
Tracer® EX2	
TREX26P	16" x 16" wheelchair, Offset Permanent Full Padded Arm, Midnight Blue Vinyl
TREX26R	16" x 16" wheelchair, Removable Desk Padded Arm, Midnight Blue Vinyl
TREX28P	18" x 16" wheelchair, Offset Permanent Full Padded Arm, Midnight Blue Vinyl
TREX28RF	18" x 16" wheelchair, Removable Full Padded Arm, Midnight Blue Vinyl
TREX28R	18" x 16" wheelchair, Removable Desk Padded Arm, Midnight Blue Vinyl
TREX20R	20" x 16" wheelchair, Removable Desk Padded Arm, Midnight Blue Vinyl
TREX28FF	18" x 16" wheelchair, Fixed Frame, Permanent Full Padded Arm, Midnight Blue Vinyl
Tracer® SX5	
TRSX56FB	16" x 16" lightweight wheelchair, Flipback Desk Padded Arm, Black Nylon
TRSX58FBF	18" x 16" lightweight wheelchair, Flipback Full Padded Arm, Black Nylon
TRSX58FB	18" x 16" lightweight wheelchair, Flipback Desk Padded Arm, Black Nylon
TRSX50FB	20" x 16" lightweight wheelchair, Flipback Desk Padded Arm, Black Nylon
	(300 lb. weight capacity)
Tracer® SX5 Recliner	
TRSX5RC6	16" x 16" Recliner, Removable Desk Padded Arm, Midnight Blue Vinyl
TRSX5RC8	18" x 16" Recliner, Removable Desk Padded Arm, Midnight Blue Vinyl
Tracer® IV	
T420RDA	20" x 18" heavy-duty wheelchair, Removable Desk Padded Arm, Midnight Blue Vinyl
Γ420RFA	20" x 18" heavy-duty wheelchair, Removable Full Padded Arm, Midnight Blue Vinyl
Γ422RDA	22" x 18" heavy-duty wheelchair, Removable Desk Padded Arm, Midnight Blue Vinyl
Γ422RFA	22" x 18" heavy-duty wheelchair, Removable Full Padded Arm, Midnight Blue Vinyl
Γ424RFA	24" x 18" heavy-duty wheelchair, Removable Full Padded Arm, Midnight Blue Vinyl
Γ424RDA	24" x 18" heavy-duty wheelchair, Removable Desk Padded Arm, Midnight Blue Viryl
T4X22RDA	22" x 18" heavy-duty wheelchair, Removable Desk Padded Arm, Midnight Blue Vinyl
	(450 lb. weight capacity)
T4X24RDA	24" x 18" heavy-duty wheelchair, Removable Desk Padded Arm, Midnight Blue Vinyl
	(450 lb. weight capacity)



Yes, you can."

**Invacare Corporation** www.invacare.com

#### USA

One Invacare Way Elyria, Ohio 44036-2125 (800) 333-6900

# Canada 570 Matheson Blvd. E. Unit 8

Mississauga, Ontario L4Z 4G4 Canada (800) 668-5324

Invacare and Tracer are registered trademarks and Yes, you can., the blue medallion design, IVC, Topaz and Jymni are trademarks of Invacare Corporation. Total One Stop Shopping is a service mark of Invacare Corporation. Ramps are manufactured by Prairie View Industries and distributed by Invacare Corporation. Velcro is a registered trademark of Velcro Industries B.V. Textilene is a registered trademark of Ludlow Corporation. Specifications are subject to change without notification.



"Invacare says 'Yes, you can",' and I couldn't have said it better myself."

© 2003 Invacare Corporation

Form No. 00-261 Rev. 10/03