

# OXYGEN CONCENTRATOR

★★★★★ PROVIDES UNLIMITED SUPPLY OF OXYGEN

## 5 LITER CONCENTRATOR

USER FRIENDLY DESIGN WITH LOW NOISE LEVEL

### SPECIFICATIONS

- Ultra-Quiet and lightweight design
- Low noise level: 46dB(A)
- Oxygen Flow: 5 LPM
- Flow Range: 0.5-5LPM
- Oxygen Concentration: 0.5-5 LPM, 93% +/-3%
- Maximum Output Pressure: 5.5 PSI +/- 0.5 PSI
- Pressure Release Mechanism: 36 PSI +/- 7 PSI
- MCU intelligent control
- Long-life compressor with 3 year warranty and maintenance-free use



5 Liter

#### PRODUCT WEIGHT

35 lbs.

#### PRODUCT DIMENSIONS

13" X 10" X 21"



### SAFETY FEATURES



Current overload or line surge shutdown



High temperature compressor shutdown



High pressure alarm shutdown



Low pressure alarm shutdown



Low oxygen concentration alarm

MODEL # OC5L