


FEATURES

- Separates into two pieces for easy carrying
- Rotation of pin either unlocks or locks the ramp sides together without extra tools
- Folds, carries like a suitcase
- Patented hinge eliminates pinch point
- Closure strap locks ramp panels together
- Lightweight, easy-to-handle & set up
- Durable welded construction
- Safety guide indicates if ramp is on a safe slope
- Full platform provides excellent stability and accommodates wheelchairs and scooters with various wheel configurations
- Anti-slip, high traction surface
- Safety DVD and steel security pins included
- Fits in the trunk of most cars for easy transportation



Item Number	Length	Width	Folded Size	Weight	Capacity*
WCR530	5'	30"	32-1/4" x 16" x 9"	34 lbs.	800 lbs.
WCR630	6'	30"	38-1/4" x 16" x 9"	40 lbs.	800 lbs.
WCR730	7'	30"	44-1/4" x 16" x 9"	44 lbs.	800 lbs.
WCR830	8'	30"	50-3/4" x 16" x 9 1/2"	52 lbs.	800 lbs.
BAR536	5'	36"	32-1/4" x 17 1/2" x 11 1/2"	50 lbs.	800 lbs.
BAR636	6'	36"	38-1/4" x 17 1/2" x 11 1/2"	58 lbs.	800 lbs.
BAR736	7'	36"	44-1/4" x 17 1/2" x 11 1/2"	67 lbs.	800 lbs.
BAR836	8'	36"	50-3/4" x 17 1/2" x 11 1/2"	76 lbs.	800 lbs.

Not sure which ramp length you need? Refer to pages 13-14 to help you.

* Maximum ramp capacity: 400 lbs. (one axle) – 800 lbs. (two axles).