

This Universal Cover Creates a Raised, Defined Perimeter to Promote a Safe Bed Environment, Without the Need for Confining Side Rails.

NEW

1P
ACK



COMPARE TO
Compare to Posey® 5750 and 5750NC

Defined Perimeter Mattress Cover

Item #	Product Dimensions	UOM
14333	34" - 39" (W) x 80" (L) x 6" - 8" (H)	1/ea
14333-42	39" - 43" (W) x 80" (L) x 6" - 10" (H)	1/ea

- The Defined Perimeter Mattress Cover creates a raised rail, defined perimeter for enhanced fall prevention, without using patient restraints.
- Open in mid-section for safe and easy patient egress and ingress.
- Built-in medium density foam bolsters safely and comfortably define edges of bed for patient and prevent falling.
- Easy-to-install and secure.
- Breathable Nylon makes this a truly universal cover, compatible with both foam and powered air surfaces.

specifications	
Fits Mattresses	39-43" (W) X 80" (L) X 6"-10" (H)
Foam Sections	5½" (W) x 29½" (L) x 4"(H)
Shipping Weight	5 lbs.
Cover Material	Nylon
Foam Density	Medium
Warranty	Limited 1 Year