


# RETRACTABED® Series

COMFORT WIDE™ (Expandable)

PUSH-BUTTON SLIDE TECHNOLOGY  
WITH LOCKING POSITIONS 35-39-42  
INCHES WIDE



EVAC POSITION™

PASS THROUGH DOORS NO MATTER  
THE ASSIST DEVICE ATTACHED

RETRACTABED®  
Model

- **Clinical-Contour®** – exclusive feature that reduces skin shear and the repositioning of residents
- **Evac Position™** – extra position to make the bed thinner so if a side bar is attached it can still pass through doors less than 42 inches wide (**expandable only**)
- **Comfort Wide™** – bed expands to 42 inches wide without the use of tools for more comfort and flexibility (**expandable only**)
- **Staff Lockout Control** – standard feature built into pendant to prevent unnecessary operation
- **Wall Bumpers** – soft non-marking bumpers come standard to protect wall from damage
- **Retracting Bed Frame** – reduces potential falls and increases staff productivity
- **Embedded Bed Controls** – optional fixed controls in assist bar or head rail reduces replacement handsets and increases staff productivity
- **SureLock** – optional one step lock for improved ease of use while locking casters

## RetractaBed® Model SPECIFICATIONS

Deck Height Low .....	9"
Deck Height High .....	25"
Expandable Sleep Surface.....	35"w, 39"w, 42"w x 80"l
Fixed Sleep Surface .....	35"w x 80"l
Maximum Safe Working Capacity .....	500 lbs
Product Weight .....	225 lbs

### LIMITED WARRANTY

15 years on bed frame mechanics	1 year on wood, casters & handsets
5 years on drive system	Lifetime warranty on all welds



**FALL RISK**  
COMPETITIVE  
BRAND

**FALL PREVENTION**  
MED-MIZER  
RETRACTABED

## ENDBOARD OPTIONS



## ACCESSORIES



PROUDLY MADE IN THE U.S.A.

(877) 867-7365  
**MED-MIZER.COM**

Approvals: IEC60601-1:2005 2nd edition, ANSI/AAMI ES60601-1:2005 2nd edition, and CAN/CSA-22.2 No 60601-1:2008



is at the forefront of innovating medical technology for post acute care. We are focused on partnering with you to help address healthcare challenges in the following key areas—




**FALLS PREVENTION**




**BARIATRIC CARE**



**MOBILITY**



**REHABILITATION**



**WOUND CARE**