

# RetractaBED™ Series

---



## RetractaBed Fixed MODEL

- **Retracting Bed Frame** — reduce potential falls and improve staff productivity
- **Clinical-Contour®** — exclusive feature that reduces skin shear and the repositioning of residents
- **Bolt on Width Kit** — bed can expand to 42 inches wide with a bolt on kit for more comfort and flexibility
- **Retracting Bed Frame** — reduces potential falls and increases staff productivity
- **Wall Bumpers** — soft non-marking bumpers come standard to protect wall from damage
- **SureLock** — optional one step lock for improved ease of use while locking casters

**RetractaBed® Fixed Model**  
**SPECIFICATIONS**

Deck Height Low .....	9"
Deck Height High .....	25"
Sleeping Surface .....	35"w x 80"l
Maximum Safe Working Capacity .....	500 lbs
Product Weight .....	225 lbs

**LIMITED WARRANTY**


15 years on bed frame mechanics	1 year on wood, casters & handsets
5 years on drive system	Lifetime warranty on all welds



**ENDBOARD OPTIONS**



**ACCESSORIES**



PROUDLY MADE IN THE U.S.A.

(877) 867-7365  
**MED-MIZER.COM**

Approvals: IEC60601-1:2005 2nd edition, ANSI/AAMI ES60601-1:2005 2nd edition, and CAN/CSA-22.2 No 60601-1:2008



is at the forefront of innovating medical technology for post acute care. We are focused on partnering with you to help address healthcare challenges in the following key areas—

