



**Smilez**  
PEDIATRICS

Made for Kids!

Loved by Families!



# PEDIATRIC STROLLER

305-363-5770

[www.smilezpediatrics.com](http://www.smilezpediatrics.com)

**The SMILEZ Adaptive Pediatric Stroller** is a versatile, long-term mobility solution designed to provide dependable positioning, comfort, and ease of use across home, school, and clinical settings. Built on a durable umbrella-style frame, it offers excellent maneuverability while supporting children with varying mobility needs.

### Key Highlights

- Adjustable seating angles for proper positioning and comfort
- Umbrella-style frame for easy handling and transport
- Large rear wheels for improved stability and maneuverability
- Multiple seat sizes with high weight-capacity options
- Suitable for families, schools, and clinical environments
- Optional accessories available for customization



### Colors to chose



## The Compact Commuter



Designed with tight spaces in mind, the **SMILEZ** Adaptive Pediatric Stroller's intuitive one-handed fold transforms it into a compact, self-standing package in one smooth motion. This hassle-free storage is ideal for families frequenting busy cafes, riding on trains, or managing limited luggage space in a hotel room.

## The Smilez Difference

### Color Options Available in Every Size

Unlike most special needs strollers that limit color choices, SMILEZ offers multiple color options across all stroller sizes—so families don't have to sacrifice style for function.

### Growth-Ready Design for Long-Term Use

SMILEZ is built to grow with your child, supporting longer use over time—so families aren't forced to replace the stroller as quickly as they would with standard models.

### Built-In Support Harness for Safer, More Secure Seating

Equipped with a supportive harness system that helps keep your child safely positioned and secure—something most everyday strollers and many standard medical strollers don't provide.

### Free Accessories Included (No Add-On Cost)

SMILEZ includes essential accessories at no extra charge—saving families money compared to other strollers, where accessories are often sold separately.

## Add-Ons and Accessories



Designed to grow with the child, the **SMILEZ™** stroller is available in multiple seat sizes with broad weight capacities, making it suitable for a wide range of pediatric and adolescent users. Color options are available across all sizes, and optional accessories including trays, head support, storage, and sun protection allow facilities and families to customize the stroller to their specific needs.

## Specifications

### SEAT SIZES

	12 inches	14 inches	16 inches	18 inches
Back Height	24"	24"	24"	24"
Seat Depth	9"-13"	10"-14"	11"-15"	12"-16"
Seat-to-Back Angle	80°- 95°	80°- 95°	80°- 95°	80°- 95°
Seat Angle	20°	20°	20°	20°
Seat Width	12"	14"	16"	18"
Seat-to-Footplate Height	4"-20"	4"-20"	4"-20"	4"-20"
Seat-to-Floor Height at 90° (rear/front)	19"/22"	19"/22"	19"/22"	19"/22"
Overall Width (outside rear wheels)	20.5"	22.5"	24.5"	26.5"
Product Weight	33 lb	34 lb	35 lb	36 lb
Weight Capacity	75 lb	100 lb	170 lb	250 lb
Weight Capacity for Transit	75 lb	100 lb	165 lb	165 lb
Front Wheels	8"	8"	8"	8"
Rear Wheels	10"	10"	10"	10"