







PEDIATRIC STROLLER

The SMILEZ™ Adaptive Pediatric Stroller is a durable, umbrella-style mobility stroller designed to provide reliable positioning, comfort, and maneuverability in home, school, and clinical environments. An adjustable seat-to-back angle (80°-95°) allows for customized positioning, while padded, removable upholstery and a secure five-point harness ensure safe and supported seating.



Specifications

	12 inches	14 inches	16 inches	18 inches
Back Height	24"	24"	24"	24"
Seat Depth	9"-13"	10"-14"	11"-15"	12"-16"
Seat-to-Back Angle	80° - 95°	80° - 95°	80° - 95°	80° - 95°
Seat Angle	20°	20°	20°	20°
Seat Width	12"	14"	16"	18"
Seat-to-Footplate Height	4"-20"	4"-20"	4"-20"	4"-20"
Seat-to-Floor Height at 90° (rear/front)	19"/22"	19"/22"	19"/22"	19"/22"
Overall Width (outside rear wheels)	20.5"	22.5"	24.5"	26.5"
Product Weight	33 lb	34 lb	35 lb	36 lb
Weight Capacity	75 lb	100 lb	170 lb	250 lb
Weight Capacity for Transit	75 lb	100 lb	165 lb	165 lb
Front Wheels	8"	8"	8"	8"
Rear Wheels	10"	10"	10"	10"
Color options				

Features

Large 10" rear wheel enhancement stability and ease of movement, with wheel locks for added safety during transfers. Available in multiple seat sizes with broad weight capacities, the SMILEZ™ stroller is built to grow with the child.