



Pediatric Wheelchair


Made for Kids! Loved by Families!
305-363-5770 www.smilezpediatric.com

The *SMILEZ Pediatric Wheelchair* is a lightweight, foldable mobility solution designed specifically for pediatric users and built to meet the needs of dealers, therapists, and facilities. Constructed with a durable steel frame and robust nylon upholstery, it delivers dependable performance while remaining easy to transport and store.

Accessories Included



Cushion & Seat Belt



Anti-Tippers



Smilez on the Move

The Smilez Kids Wheelchair sets a new standard in pediatric mobility. Through its durable steel frame construction and highly adjustable seating systems, it addresses the complex physiological needs of growing children while ensuring long-term durability and ease of use for caregivers

The Smilez Difference

Color Options Available in Every Size

Unlike many special needs wheelchairs that offer limited or no color variety, **SMILEZ wheelchairs come in multiple color options across all sizes**, allowing families to choose a look their child actually loves—without sacrificing function or support.

Built-In Support Buckle for Safer, More Secure Seating

SMILEZ wheelchairs include a **built-in support buckle / secure seating system** that helps keep the user safely positioned and stable. This added security supports children who need extra safety and posture control—something most basic wheelchairs do not include by default.

Free Accessories Included (No Add-On Cost)

Unlike many wheelchairs where key add-ons cost extra, **SMILEZ includes essential accessories for free**, helping families save money while getting everything they need from day one—no surprise upgrades or extra purchases required.

Adjustable Footrests and Armrests

SMILEZ wheelchairs feature adjustable footrests and armrests, allowing for better fit, comfort, and support as the child grows. These adjustments improve positioning, reduce discomfort, and create a more customized ride compared to standard chairs with fixed components.

Features



Height Adjustable Push Handles



Padded flip-back armrests



Built-in bus transit tie-downs



Height-adjustable elevating legrests



Specifications >>

Multiple kid-friendly color options



Seat Width

12"

14"

Seat Depth

12"

16"

Back Height

14"

16"

Seat Height Front

18"

18.5"

Seat Height Back

17"

17.5"

Overall Open Width

18.5"

21"

Overall Length with Riggings

36"

41"

Weight w/o Riggings

28 lb

36.5 lb

Weight Capacity

150 lb

150 lb