



## PEDIATRIC WHEELCHAIR

The SMILEZ™ Pediatric Wheelchair is a lightweight, foldable mobility solution designed specifically for pediatric users and built to meet the needs of dealers, therapists, and facilities. Constructed with a durable steel frame and robust nylon upholstery, it delivers dependable performance while remaining easy to transport and store.



### Specifications

	12 inches	14 inches
Seat Width	12"	14"
Seat Depth	12"	16"
Back Height	14"	16"
Seat Height Front	18"	18.5"
Seat Height Back	17"	17.5"
Overall Open Width	18.5"	21"
Overall Length with Riggings	36"	41"
Weight w/o Riggings	28 lb	36.5 lb
Weight Capacity	150 lb	150 lb

### Features

- Easily foldable lightweight design
- Cushion, seat belt, and anti-tippers included
- Padded flip-back armrests
- Height-adjustable elevating legrests
- Built-in bus transit tie-downs
- Multiple kid-friendly color options

