

# The Safety Sleeper (TSS)



## MODEL COMPARISON



Made in the USA



Machine washable & dryable



Fire retardant, medical grade fabric



Easy assembly, no tools required



Less than 50 pounds  
\*Excluding Queen Model\*



The most versatile bed on the market, for home & away



Safety lab tested & certified



Compatible with most beds or bed frames

### 200 MODEL



- ✓ Set of 4 frame pads
- ✓ 2 stabilizing straps
- ✓ Air mattress for travel
- ✓ Wheeled suitcase
- ✓ 8" memory foam waterproof mattress
- ✓ Freestanding frame

- ✓ 2 Coverlet Sheets
- ✓ One side entry

\*Customizable Model

### 300 MODEL



- ✓ Set of 4 frame pads
- ✓ 2 stabilizing straps
- ✓ Air mattress for travel
- ✓ Wheeled suitcase
- ✓ 8" memory foam waterproof mattress
- ✓ Freestanding frame

- ✓ Side entry
- ✓ Access points (1 each side)
- ✓ 2 Waterproof Coverlet Sheets
- ✓ End of bed entry with removable vinyl and mesh panels

### 400 MODEL



- ✓ Set of 4 frame pads
- ✓ 2 stabilizing straps
- ✓ Air mattress for travel
- ✓ Wheeled suitcase
- ✓ 8" memory foam waterproof mattress
- ✓ Freestanding frame

- ✓ Convertible entry
- ✓ Lift Compatible
- ✓ Access points (1 each side)
- ✓ 2 Waterproof Coverlet Sheets
- ✓ End of bed entry with removable vinyl and mesh panels

#### Choice of 6 colors:



Heather Grey



Hot Pink



Navy Blue



Royal Purple



Sail Blue



Shimmering Sand\*

\*shown in photos above



# MODEL OVERVIEW



#### FRAME PADS & BUILT IN PADDING:

Dense foam, encased in fabric, wraps around the top and sides of the metal frame to protect against seizures and self-harming behaviors. Easy zip on/off for removal. Machine washable.



#### COVERLET SHEETS:

Replace the need for a fitted sheet by zipping into TSS to protect against entrapment zones. Coverlets are only accessible from the end of the bed zipper and can be easily removed to launder separately from the entire enclosure.



#### STABILIZING STRAPS:

Straps are attached to TSS frame and encompass your traditional bed frame to prevent tipping or moving. Recommended to counteract rambunctious activity.



#### MATTRESS:

TSS comes standard with an 8" waterproof memory foam mattress.



#### REFLECTIVE PULL-TAB VELCRO SECURING SYSTEM:

For ease of opening and safeguarding, all doors are secured with a reflective handle. All zippers extend over webbing to prevent any zipper head access from inside the bed.



#### ACCESS POINTS:

One 1-1/8 inch metal grommet installed on both side entries (2 total) of the enclosure to serve as access points for use with medical tubing. Alternate access points available upon request.



#### WHEELED SUITCASE & AIR MATTRESS:

TSS comes packaged in its own hard-shell, wheeled suitcase. Included for ease of travel is an air mattress with a built-in electric pump.



#### END OF BED ENTRY WITH REMOVABLE VINYL AND MESH PANELS:

This comes standard with 300 and 400 models. Optional upgrade with 200 model.

	Overall Dimensions	Headroom	Included Mattress
<b>Twin</b>	74.5"L x 41.5"W x 44"H	36"	74.5"Lx37"W
<b>Twin XL</b>	79.5"L x 41.5"W x 44"H	36"	79.5"Lx37"W
<b>Full</b>	74.5"L x 56.5"W x 44"H	36"	74.5"Lx52"W
<b>Queen</b>	79.5"L x 62.5"W x 44"H	36"	80"Lx58"W

**Customizable Options:** \*200 model

- Waterproof coverlet sheets
- Inside pocket
- End of bed entry with removeable vinyl
- Convertible entry
- Access points

Patented design: US Patent #7,971,596 B2 Registration #3011389091 - Owner/Operator #10048434.

