


INVACARE®  
**9000 XDT**  
WHEELCHAIR



The **Invacare® 9000 XDT Wheelchair** is designed for individuals who demand an extra wide, heavy-duty or tall chair with seat-to-floor heights ranging from hemi through tall. The high-performance 9000 XDT wheelchair provides a 350 lb. weight capacity and the components are interchangeable with the other 9000 series wheelchairs.

## Features

- Durable, low-maintenance carbon steel frame is long-lasting
- Tall frames are designed to accommodate taller individuals and provide a more comfortable ride
- Multiple frame color variations for customizing
- Triple-axle positions allow wide variation of seat-to-floor heights
- Urethane rear tires mounted on "no flex" wheels offer superior performance
- Urethane casters for better ride and lighter weight
- Low-maintenance, dual-sealed precision bearings
- Heavy-duty inner liners keep seat and back from stretching

HCPCS Code: K0007



Invacare Corporation

1. Triple-axle positions
2. Urethane casters
3. Bariatric sizes available
4. Seating available
5. Builder options

Model #	Overall Width	Overall Height	Overall Depth w/o riggings	Folded Width	Seat Dimensions	Seat-to-Floor Height (2" increments)	Back Height
9XDT	23"-29"	36"	29"	13"	16"-20", 22" x 16"-20"	17.5"-21.5"	14"-19"
	Caster Height	Wheel Height	Weight w/o Rigging	Weight Capacity	Shipping Weight	Carton Dimensions	Limited Warranty
9XDT	8"	24"	37 lb.	350 lb.	55 lb.	34"x13"x 39"	5 Years

©2014 Invacare Corporation. All rights reserved. Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Specifications are subject to change without notification. Form No. 11-229 Rev. 04/14 131071



All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further detail.