



Redefining standards in care



A new adjustable electric bed

Safe, Versatile, Efficient



Yes, you can.®

Technical data ▼ Blue: Indicates Item Number



Bed Deck Length ▼



Bed Deck Width ▼



Overall Length (incl. head/foot boards) ▼



Maximum Safe Working Load ▼



Product Weight ▼



Minimum/Maximum Bed Deck Height ▼

CS600-36	80"	36"	88"	600 lbs	217 lbs	7" - 30"
CS600-EXP	80"	36", 39", 42"	88"	600 lbs	243 lbs	7" - 30"




Maximum Back Angle ▼



Maximum Knee Angle ▼



Maximum Reverse Trendelenburg ▼



Length Extension ▼



Width Extension ▼



ACP (Attendant Control Panel) ▼

CS600-36	60°± 5°	30°± 5°	17°± 5°	80" - 84" w/accessory	Optional w/accessory	Optional
CS600-EXP	60°± 5°	30°± 5°	17°± 5°	80" - 84" w/accessory	Integrated	Optional

Bed end finishes



Solar Oak

36"
CS600BE36SO
42"
CS600BE42SO



African Walnut

36"
CS600BE36AF
42"
CS600BE42AF



Biltmore Cherry

36"
CS600BE36BC
42"
CS600BE42BC

Frame finish



Light Gray

Bed end style



Amherst

Please note: colors may vary slightly from those displayed above.

Quarter length rail options

Standard width, CS600-36

Head section, Qty: 2
CS600HALFRAILH

Foot section, Qty: 2
CS600HALFRAILF

Expandable width, CS600-EXP

Head section, Qty: 2
CS600HALFRAILHEXP

Foot section, Qty: 2
CS600HALFRAILFEXP

Assist rail options

Standard width, CS600-36

Assist rail, Qty: 1
CS600ASTRAIL36

Expandable width, CS600-EXP

Assist rail, Qty: 1
CS600ASTRAILEXP

Other accessories

Length extension
CS600LENGTHEXP

Trolley
CS600TROLLEY

Battery backup
CS600BATBCK

Trapeze
CS600TRAPEZE

Width extension
CS600WIDTHEXP

Underbed light
CS600UNDLIGHT

ACP (Attendant Control Panel)
CS600ACPPAD

Tel: (800) 333-6900
www.invacare.com

Invacare Corporation
Worldwide Headquarters
1 Invacare Way
Elyria, OH 44035-4190



Yes, you can.®

©2024 Invacare Corporation. All rights reserved. Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Specifications are subject to change without notification. Please refer to the User Manual for instructions for use. Form: 23-037 Job: 240136US-04/2024