



Advantages

Adjustable Arm & Foot Rest
The flip-arm and foot rest allows for easier patient access and patient comfort

2 Wheels
The base of the scale has two wheels that provide easy mobility and allows for convenient placement



Model	Shipping Dims. (L x W x H) / Cube	Weight (LB)	UPC Code
445KL	30.25" x 25.5" x 45" / 20.09	81.5 lb	892076002046



Specifications

- Capacity: 440 lb / 200 kg
- Graduation: 1/4 lb / 100 g
- Seat Size: 18 1/2" (w) x 17 1/4" (d) x 15 3/4" (h) / 470 mm (w) x 438 mm (d) x 400 mm (h)
- Wheels: 2
- 2 Year Limited Warranty