

ICARE IC100 PRODUCT DETAIL



BED CODES & SIZES

TWIN XL

Code: IC100TXL

Mattress Surface Size: 80" x 35"

Base Size: 85" x 38"

KEY FEATURES

- 440lbs Safe Working Load
- Optional Headboards and Footboards
- Compatible with Icare Bed Accessories
- 5 Year Warranty

COMPATIBLE ACCESSORIES AVAILABLE

- High Side Rail
- Low Side Rail
- U-Assist Rail
- Over Bed Pole
- Headboard / Footboard
- Slide-Out Side Rail
- Over Bed Tables
- Bedside Folding Shelf
- Bedside Tables
- Blanket Cradle

TECHNICAL SPECIFICATIONS

- Weight rating – The max weight limit is 440lbs includes any model of Icare mattresses.
- Fabric – Hard wearing commercial grade fabric. Vinyl is treated with antibacterial properties.
- Mattress Surface – Strong metal mesh under mattress allowing increased airflow.
- Wheels – large lockable heavy duty swivel castors on each corner.
- Fixed base height – 14" from floor to mattress platform.
- Under Bed clearance – 7"