





Client First

Headquartered in High Point, North Carolina, we strive to provide the highest quality products and customer service through our team of associates and three distribution centers in the US.



Charme Inc. is incorporated in Delaware, US and operates across the US, Europe and Asia. Through starsleep™, Charme is committed to revolutionize the world of smart bed industry with its strong manufacturing capabilities and global distribution network. Our team of sleep specialists and engineers are dedicated to technological innovation and providing the best possible sleeping experiences.



Like nothing else on the market today, Orin provides care, comfort, security and independence to the highest standard through our innovative product design.



Our Customer Profiles



Diabetes

Pregnancy

Acid Reflux

Sleep Deprivation

Cancer

Stroke

Arthritis

Wheelchair Bound

All-In-One Features







4D Air Fiber Mattress Washable & Breathable



German Technology

4









Head & Foot Articulation






One Touch Flat



Sleep Report
















Heating



WATCH VIDEO







Ergonomic Design

Like nothing else on the market today, Orin provides care, comfort, security and independence to the highest standard through our innovative product design.

Orin Rotation Bed Key Features

- Multifunctional Design
- Silent & Durable German Motor
- One-Touch 90° Rotation Bed
- Head (90°) and Foot (35°) Elevation
- Air Chamber Massage
- Heating
- Air Fiber Breathable Mattress
- Sideboards & Safety Belt For Protection
- US FDA Registered For Bed & Mattress
- Sleep & Health Monitoring Apps
- 170kg Weight Capacity









Multifunctional Design One-Click Control

90° rotation helps users to get in and out of bed easily.



Rotate 90° (easy to get out)

Step 3

Steady Rotation


Orin Rotation Bed Dimension

2 Bed Models for Selection

Left Rotation | Right Rotation







Air Massage with Heating Option

The massage airbags consist of four independently controllable airbags that are installed inside the mattress for use. Three massage modes can be achieved by inflating and deflating the airbags. The waist airbags can also function as lumbar support. If there is no operation, the airbags will automatically stop running after 45 minutes.

The heating function is provided by a 29V heating pad. If left unattended, the heating pad will automatically stop after running continuously for 60 minutes.

3 Modes To Offer Personalized Care

4-Zone Massage

40°C Heating




Air Chamber Massage



Heating










German Technology

High Efficiency / Silent / Safe & Stable









Head & Foot y Articulation







> 25000order Durable

> 170kg Stable Load-Bearing Flame Retardant Medical Standard






Ensures Normal Use During Power Outage

Improve care and reliability.








Spine-Protecting High-Density Mattress

Designed to fit every body type. Breathable, silent & non-disturbing. Instantly boost sleep quality.







Ergonomic Design



Zero-G





Air Chamber

◀ High Density Foam

Allergy-Fighting Fabric

Allergy-Fighting Fabric



Application Scenarios

Caregiver Group

Medication Record

Diet Record

Chronic Disease Management

Overnight Risk Warnings





Non-intrusive Sleep Tracking Technology

Real Time Alarm & Notification For The User & Family Members








Sleep Report



Out of Bed Alarm



Orin-Accessories

Safety Protection



Step 1 Lock Safety Belt


Step 2Adjust Tightness

Orin-Accessories Safety Protection Sideboards Protection

Orin-Accessories Overbed Table

- Multi-Home Scenarios
- Mobile, Sturdy & Durable
- Practical & Space Saving
- Air-Spring Assist, Easy Operation
- 75cm to 115cm Height Range
- 30kg Weight Capacity



















Orin-Accessories Transfer Chair

- Ergonomic Design
- Electric Easy Height Adjustment
- Easy Access To Toilet
- Cleanable
- Chargeable Battery
- 25cm Height Adjustment
- Mute and Lockable Casters
- 135kg Safety Weight Load

















Orin-Accessories Pillow / Home Textiles







Breathability









Our Factory

Starsleep's products are produced leveraging Charme's global supply chain system, incorporating Italian design, German technology, Swiss innovation, and North American quality control standards. By combining materials and production cost advantages from various regions in Asia, we optimize the manufacturing process. We have always advocated the concept of a high-quality lifestyle, embracing a life philosophy of 'Freedom, Balance, Wisdom and Taste,' and embodying the brand principles of 'respect, integration, innovation, and gratitude.' Starsleep not only focuses on high-quality materials, practical products, and reasonable pricing but also strives to provide users with a more practical, innovative, and comfortable home experience, making the wisdom and enjoyment of life within easy reach.



j

2019193 2019139





⊀r starsleep