

# Supernal<sup>®</sup>

RECLINER BEDS



MAXIMIZE YOUR  
**INDEPENDENCE** WITH OUR  
**HI-LOW BEDS**

# Supernal

HI - LOW



Supernal Hi-Low Dual King

## FEATURES:

- Head (65°), Foot (35°), and Hi-Low adjustable
- Hugs the wall so you can reach your nightstand
- Super low deck height for easy self transfer (10.5")
- Excellent height range for flexibility of use (10")
- Standard weight capacity (400/500 lbs.)
- Cradles your pillow with a European style head tilt
- Backlit wireless hand control for visibility
- Variable massage settings
- Emergency power down



The **Supernal Recliner** is an optional 2 function version of this bed (deck height 11.5"). Pair it with the Hi-Low for a **Dual King** setup.

Hospital bed features without the hospital bed look.  
Available sizes: Twin 80, Full 80, Queen, and Dual King.

# Supernal 5



## FEATURES:

- Head (60°), Foot (35°), and Hi-Low Adjustable
- Trendelenburg (10°), Reverse Trendelenburg (10°), and Cardiac Chair (70°)
- Low deck height for easy self transfer (12")
- Maximum height range for greatest flexibility (12")
- Standard weight capacity (300 lbs.)
- Touch & Glow wireless hand control
- Variable massage settings



Hospital bed features without the hospital bed look.  
Available sizes: Twin 80, Full 80, Queen, and Dual King.

# Supernal<sup>®</sup>

## SLEEP SYSTEM

The Supernal Sleep System low-profile mattresses are hand selected as perfect matches for our Supernal beds—and most adjustable beds. Unique border edge protection prevents slipping.

Two mattress styles meet a wide range of needs:

- **Ascent** - our firmer mattress
- **Soft Touch** - our softer memory foam mattress



- Two mattress cover options
- Zipper construction for easy cover replacement



Natural bamboo quilted cover



Stretch fabric cover (waterproof)

### AVAILABLE OPTIONS

Customize your Supernal bed with the following options (availability depends on model):

- Half Rails
- Bamboo Rail Covers
- Headboard
- 5" Locking Caster (Hi-Low only)

